

**California Environmental Protection Agency**



# **Air Resources Board**

**PROCEDURE FOR DETERMINATION OF EQUIVALENT SOLVENT RED 26 DYE  
CONCENTRATION IN DIESEL FUELS BY PORTABLE VISIBLE SPECTROMETER**

**Standard Operating Procedure No. MLD 154**

Revision No. 1.0

Effective date: April 2012

Southern Laboratory Branch  
Monitoring and Laboratory Division  
9528 Telstar Avenue  
El Monte, CA 91731

This report has been reviewed by the staff of the California Air Resources Board and approved for publication. Approval does not signify that the contents necessarily reflect the views and policies of the Air Resources Board, nor does the mention of trade names or commercial products constitute endorsement or recommendation for use.

# **SOP No. 154 - PROCEDURE FOR DETERMINATION OF EQUIVALENT SOLVENT RED 26 DYE CONCENTRATION IN DIESEL FUELS BY PORTABLE VISIBLE SPECTROPHOTOMETER**

## **1 Introduction**

- 1.1 This document describes a method of determining the concentration of red dye levels equivalent to 0.1 to 20 mg/L of Solvent Red 26 in commercially available diesel and burner fuels using visible spectroscopy.
- 1.2 This method is based on ASTM 6756-02.

## **2 Method Summary**

The diesel samples are transferred into a cuvette, which is placed into the light path of the instrument. A beam of visible light is imaged through the sample onto a detector, and the detector response is determined. Wavelength ranges which correlate strongly with the red dye concentration are selected for analysis using bandpass filters. A multivariate mathematical model converts the absorption values to the red dye concentration.

## **3 Interferences and Limitations**

The presence of colorants other than typical diesel hydrocarbons, or the presence of red dye other than the specified types, can interfere with the accurate determination of the red dye concentration.

## **4 Instrumentation and Apparatus**

- 4.1 Petrospec DT-100C portable spectrometer
- 4.2 Sample cells (cuvettes): fused silica, glass, or polymethacrylate cells with a sample path length of 10 mm.

## **5 Reagents and Materials**

Check samples containing red dye 26 equivalent concentrations of 0.0, 5.0, and 12.0 mg/L can be prepared in-house or purchased from the instrument manufacturer.

## **6 Safety Precautions**

Standard laboratory safety procedures and equipment should be used in performing this method. For example, safety glasses and gloves should be worn. All standard and sample preparation should be done in the fume hood. Diesel fuel contains compounds known to be toxic and carcinogenic.

## **7 Samples**

