

California Environmental Protection Agency



SOP MLD 057

**STANDARD OPERATING PROCEDURE FOR THE  
DETERMINATION OF  
1,3-BUTADIENE AND BENZENE IN AMBIENT AIR BY  
CAPILLARY COLUMN GAS CHROMATOGRAPHY WITH  
PHOTOIONIZATION DETECTOR**

Engineering and Laboratory Branch  
Monitoring and Laboratory Division

First Approved Date of SOP:	January 1, 1999
Approval Date of Last SOP Amendment:	January 1, 1999
Revision Number:	1.0

DISCLAIMER: Mention of any trade name or commercial product in this Standard Operating Procedure does not constitute endorsement or recommendation of this product by the Air Resources Board. Specific brand names and instrument descriptions listed in the Standard Operating Procedure are for equipment used by the ARB laboratory.

## **Table of Contents**

1.0 SCOPE.....	1
2.0 SUMMARY OF METHOD.....	1
3.0 INTERFERENCES AND LIMITATIONS.....	2
4.0 APPARATUS .....	3
5.0 REAGENTS.....	5
6.0 INSTRUMENT CONFIGURATION AND PARAMETERS .....	6
7.0 DAILY OPERATION.....	8
8.0 DATA ANALYSIS.....	10
9.0 QUALITY CONTROL .....	10
APPENDIX I – VARIAN STAR CHROMATOGRAPHY WORKSTATION .....	33
APPENDIX II – ADDITIONAL SETPOINTS .....	48
APPENDIX III - PERKIN-ELMER NELSON TURBOCHROM DATA SYSTEM.....	49

## **List of Tables**

Table 1: Function of Valves for Varian GC 3400CX System .....	14
Table 2. List of Valve Actuators for Varian GC 3400CX System.....	14
Table 3: Pre-Injection Program Times, Relay #s, and Status .....	15
Table 4: GC Relay Program Times, Relay #s, and Status .....	16
Table 5: Estimated Limits of Detection.....	17
Table 6: Precision Measurements of 1,3-Butadiene and Benzene.....	18
Table 7: MLD057 Standards and Controls.....	19

## **List of Figures**

<b>Figure 1: Pre-Injection Program - Idle State .....</b>	<b>20</b>
<b>Figure 2: Pre-Injection Program - Purge Sample Line.....</b>	<b>21</b>
<b>Figure 3: Pre-Injection Program - Sample Loaded .....</b>	<b>22</b>
<b>Figure 4: Pre-Injection Program - Sample Loading and Cryotrap Purge.....</b>	<b>23</b>
<b>Figure 5: Pre-Injection Program - Injector Isolated for Heating .....</b>	<b>24</b>
<b>Figure 6: Post-Injection Program - Varian GC Control Started .....</b>	<b>25</b>
<b>Figure 7: Post-Injection Program - Transition and Cryotrap Pre-Heated .....</b>	<b>26</b>
<b>Figure 8: Post-Injection Program - Sample Injected onto the Column .....</b>	<b>27</b>
<b>Figure 9: Injector Temperature Conversion .....</b>	<b>28</b>
<b>Figure 10: Chromatograms.....</b>	<b>29</b>
<b>Figure 11: Multipoint Analysis for 1,3-Butadiene.....</b>	<b>30</b>
<b>Figure 12: Multipoint Analysis for Benzene .....</b>	<b>31</b>
<b>Figure 13: Control Charts .....</b>	<b>32</b>



# **SOP MLD 057**

## **STANDARD OPERATING PROCEDURE FOR THE DETERMINATION OF 1,3-BUTADIENE AND BENZENE IN AMBIENT AIR BY CAPILLARY COLUMN GAS CHROMATOGRAPHY WITH PHOTOIONIZATION DETECTOR**

### **1.0 SCOPE**

This document describes the procedures followed by Monitoring and Laboratory Division (MLD) staff to analyze 1,3-butadiene and benzene by Gas Chromatography (GC) in ambient air samples collected from the California Toxic Monitoring Network. The method was developed by staff of the Engineering and Laboratory Branch (ELB), Organics Laboratory Section (OLS). This Standard Operating Procedure (SOP), MLD057 Revision 1.0, supersedes SOP MLD051 Revision 3.0 for determining 1,3-butadiene and benzene concentrations in ambient air. The new SOP includes a multi-point calibration procedure for quantitation of butadiene and benzene. Minor changes were also made in the GC relay program, which affects the post-injection system configuration.

### **2.0 SUMMARY OF METHOD**

Ambient air is collected in a SUMMA polished stainless steel canister using a XonTech 910A sampler. The sampling procedure for Toxics samples is detailed in the Air Resources Board Quality Assurance Manual, Volume II, Appendix Q. All the operational procedures and sampling conditions for each sample are documented in the field. A record of this information is sent back to the OLS along with the sample. Upon receipt, the sample canister pressure is measured with a calibrated external pressure gauge. This information and particulars of the collection are documented in the laboratory. The sample is then analyzed according to the SOP in the laboratory.

An ambient air sample is introduced into the analytical system from a pressurized canister through stainless steel or Teflon tubing with the aid of a mass flow controller (MFC) and a vacuum system. A digital readout attached to the MFC provides a visual indication of the proper sample flow during sampling. Automated sampling of up to 16 canisters can be accomplished using the system's multi-position stream selector valve.

The sample passes through a Nafion™ dryer to remove moisture from the gas stream. It is trapped on a cryotrap at -150 degrees centigrade (°C). At this temperature, the desired components are solidified, while fixed gases, such as nitrogen (N<sub>2</sub>), oxygen (O<sub>2</sub>), and carbon dioxide (CO<sub>2</sub>), and methane (CH<sub>4</sub>) pass

through the cryotrap to the vent. The system is purged with ultrapure N<sub>2</sub> to flush sample remaining in the tubing or valving on to the cryotrap, and to remove any excess light impurities. The cryotrap is isolated and rapidly heated to 125°C, followed by injection of the sample onto a DB-VRX capillary column. After a short hold at 125°C, the trap temperature is raised to a final temperature of 190°C.

The sample mixture is separated into individual components by their interaction with the capillary column stationary phase, using temperature-programmed gas chromatography. The components eluting from the column are detected by a Photoionization Detector (PID). The target analytes, 1,3-butadiene and benzene, are subsequently identified and quantified. Routine identification of a component in a sample as 1,3-butadiene or benzene is based upon the relative retention time. The initial positive identification of a peak under the method conditions requires the use of a gas chromatograph/mass spectrometer (GC/MS) system.

### **3.0 INTERFERENCES AND LIMITATIONS**

- 3.1 Although preliminary studies have shown that 1,3-butadiene and benzene can be considered stable in stainless steel canisters, every effort must be made to analyze the sample as soon as possible. Extreme care must be taken to prevent contamination during sample collection, transportation and subsequent analysis.
- 3.2 Since the retention time of an analyte is one of the keys for identifying a component in gas chromatography, the retention times of the GC system must be reproducible.
- 3.3 Compounds with similar GC retention times may not be completely separated or may co-elute. This can cause misidentification or inaccurate quantitation.
- 3.4 Although the analytical system separates and detects other saturated and unsaturated hydrocarbons, only 1,3-butadiene and benzene are addressed by this procedure.
- 3.5 No more than 10 samples should be run consecutively without system re-calibration due to the continuous reduction in PID lamp signal.
  - 3.5.1 The loss of PID lamp signal and increased background noise can be caused by the deposition of high boiling compounds on the PID lamp lens. To reduce sample deposition, the detector temperature can be increased, but it should not exceed the maximum recommended operating temperature of 280°C. Operating at temperatures greater

than 250°C reduces the life of the internal Vespel seals. It can also **increase** background noise.

3.5.2 When the signal decreases substantially, the lamp should be cleaned or replaced. The PID lamp may be cleaned with PID lamp polishing compound followed by rinsing with water or methanol.

3.5.3 References:

“Model 703 Photoionization Detector (PID)”, Operations Manual, 116439 M, Revised March 1992, Tremetrics

“Photoionization Detector”, 03-9141118-00, July 1986, Varian

3.6 Daily baseline shifting, or the appearance of broad, extraneous “ghost” peaks, may be caused by high boiling compounds being trapped on the column. The column should be baked out prior to each set of analytical runs to remove these contaminants. The bake out temperature should not exceed the column’s maximum operating temperature of 260 °C.

3.6.1 Reference:

“1996/1997 Catalog and Technical Reference”, J & W Scientific, Inc.

3.7 Vinyl chloride is monitored by this method but is not reported due to the low ambient concentrations.

## 4.0 APPARATUS

4.1 A Varian Model 3400CX gas chromatographic system, with:

4.1.1 An automated sampler, consisting of an 8134 multi-position Stream Selector Valve (SSV) and a Mass Flow Controller (MFC) with a Control/Digital Readout module.

4.1.1.1 The MFC is mounted downstream of the SSV and the cryotrap to eliminate any contamination and to reduce dead volume in lines from sample trap.

4.1.1.2 The MFC is typically rated at 100 cm<sup>3</sup>/min at 100% full scale. The flow rate is set as a percentage of full scale. For example, a flow rate of 50 cm<sup>3</sup>/min corresponds to a setting of 50% full scale.

- 4.1.1.3 The Control/Digital Readout module is set to the side or on top of the GC.
    - 4.1.1.4 A rotometer is mounted on the GC, between the MFC and the vacuum source, to allow visual confirmation of flow.
  - 4.1.2 A Cryogenic Pre-Concentrator, containing a 150 microliter (μL) empty nickel tubing cryotrap.
  - 4.1.3 A Photoionization Detector (PID), using a 10.2/10.6 eV lamp.
  - 4.1.4 Digital flow controllers and pressure regulators installed in a heated pneumatics compartment.
    - 4.1.4.1 A pressure regulator set in parallel with the carrier gas flow controller to minimize flow upsets during valve switching.
    - 4.1.4.2 The digital flow controllers are calibrated to deliver helium (He) flows from zero to 100 cm<sup>3</sup>/min,  $\pm$  3%, with an inlet pressure of 80 psi.
  - 4.1.5 Analog pressure gauges and digital Electronic Pressure Readouts (EPR) for use in gas monitoring and diagnosing problems with the flow system.
- 4.2 A canister sampling manifold for connecting canisters to the automated sampler, using appropriate tubing and fittings.
  - 4.2.1 Examples of tubing size and material are 1/8 inch Teflon tubing, 1/16 inch stainless steel tubing, 1/16 inch nickel tubing, or 1/16 inch glass lined stainless steel tubing.
  - 4.2.2 Canisters are connected to the manifold, and the manifold is connected to the automated sampler's SSV.
- 4.3 A J&W DB-VRX 75 m by 0.45 mm i.d., with 2.55 μm film thickness, fused silica capillary column.
  - 4.3.1 Reference:

"1996/1997 Catalog and Technical Reference", J & W Scientific, Inc.
- 4.4 A continuous, self-regenerating, in-line Nafion™ sample dryer, from Perma Pure Inc.



- 4.5 A Varian GC Star Workstation, which includes an Intel compatible PC, up to four (4) Varian serial communication cards, and Varian Star Chromatography software.
  - 4.5.1 The Workstation is used for GC system configuration, sample file lists, sequence lists, and method building.
  - 4.5.2 The serial communication cards provide digital communication with the GC.
  - 4.5.3 Reference:

“Varian GC Star Workstation Manual” by Randall Bramston-Cook.
- 4.6 A Perkin-Elmer (PE) Data Station, including an Intel based PC and PE-Nelson 2700/Turbochrom™ Chromatography software, and one or more PE-Nelson Model 970 Analog/Digital Converters.
  - 4.6.1 The Data Station is used for storage of raw data files and the subsequent processing of the raw data to produce qualitative/quantitative data.
  - 4.6.2 The PE-Nelson Model 970 Analog/Digital Converter converts the GC’s analog detector signal to digital data and stores the digital raw data until it is transferred to the PE Data Station.
- 4.7 Stainless steel SUMMA passivated canisters for sample collection and standard preparation.

## 5.0 REAGENTS

- 5.1 A system blank, consisting of zero air, ultrapure air or Grade 5 ultrapure N<sub>2</sub>, in a SUMMA passivated canister that has been humidified with 150 µl of HPLC grade water.
- 5.2 A certified National Institute of Standards (NIST) hydrocarbon standard mixture containing all analytes of interest. This standard should be slightly higher in concentration than the typical sample and must be within the linear dynamic range of the GC system. [Table 7](#) lists the NIST Standards associated with the analysis of 1,3-butadiene and benzene.
- 5.3 A control sample containing 1,3-butadiene and benzene at concentrations within the calibration range of the GC System. [Table 7](#) lists the Control Standards associated with the analysis of 1,3-butadiene and benzene.

- 5.4 One high pressure gas cylinder of Grade 5 or better He for use as the GC column carrier gas.
- 5.5 Three high pressure gas cylinders of Grade 5 or better N<sub>2</sub>.
  - 5.5.1 One cylinder supplies gas to operate the pneumatic valve actuators and provides clean blank gas to port number one on the SSV. A separate gas cylinder is recommended for the valve actuators to minimize spikes on the PID from pressure pulses when valves are turned. He can be used for this purpose if necessary.
  - 5.5.2 One cylinder supplies dry gas to the Nafion™ dryer.
  - 5.5.3 One cylinder provides makeup gas to the PID.
- 5.6 A dewar of liquid nitrogen (LN<sub>2</sub>) for cooling the cryotrap and the GC column oven.

## 6.0 INSTRUMENT CONFIGURATION AND PARAMETERS

- 6.1 The analytical system's gas flow and automation configurations are shown in [Figure 1](#) through [Figure 8](#). A complete listing of the current Varian Star Workstation method, which includes all of the GC setpoints controlled by the Workstation, is given in [Appendix I](#).
- 6.2 GC Injector [Cryotrap]
  - 6.2.1 In this system, the GC injector is the cryotrap. The cryotrap utilizes a low trap volume ( $\leq 150$   $\mu$ L) and a rapid heating/cooling cycle (-150°C to 125°C in less than 0.80 minutes) for efficient transfer of sample to the capillary column.
  - 6.2.2 The set temperature for the GC Injector does **not** correspond to the actual temperature of the injector. This results from a hardware limitation in the Varian GC's on-board alphanumeric display. The actual temperature is obtained by using the Injector Temperature Conversion chart given in [Figure 9](#). The actual temperature values, as given in [Appendix I](#), are from that chart.
  - 6.2.3 GC Auxiliary

In this system, the GC Auxiliary oven is used to maintain the three Gas Sampling Valves (GSVs) at a constant temperature.

### 6.3 GC Column Oven

- 6.3.1 This oven contains the chromatographic column. Reproducible temperatures in the column oven under isothermal and programmed conditions are very important to producing reproducible retention times for components. This is critical for the qualitative aspects of the method.
- 6.3.2 The oven temperature conditions are directly involved with the interaction between all components and the stationary phase of the column. The temperature program can change the resolution between components, which affects the qualitative and quantitative aspects of the method.
- 6.3.3 The [GC Column A Parameters](#) and the [GC Column B Parameters](#) in the method DO NOT alter or control any setpoints in the GC. Listing a column as installed allows the pressure readings from the digital Electronic Pressure Readouts (EPR) to be displayed on both the GC's on-board alphanumeric display and the Star screen of the associated Workstation.

### 6.4 GC Detector(s)

- 6.4.1 The instrument has both a PID and an ECD installed in a series configuration. The ECD is not used for this analysis.
- 6.4.2 The [Attenuation](#) parameter in the method has no effect on the signal output to either the Star or PE-Nelson Chromatography Data Acquisition systems.

### 6.5 GC Autosampler

This section of the method designates the 8134 multi-position Stream Selector Valve as the input source for samples into the system.

### 6.6 GC Relays [Post Injection Relay Time Program]

This portion of the method contains the time setpoints for activating relays, which in turn control the various valves used in the GC. These actions take place after the Pre-Injection Program (PIP), which is described [below](#).

## 6.7 GC Stripchart

There is no stripchart recorder attached to the instrument. Printed output is obtained through the Star GC Workstation and the Star Chromatography software.

## 6.8 GC Gas Flows

The current gas flows and regulator pressures are shown in .

## 6.9 Mass Flow Controller

The sample volume for the column injection is fully automated by the use of 8134 SSV Relay Time Program that is part of the Varian GC Star Workstation software. The function of the valves in the Varian 3400CX GC system are shown in [Table 1](#). The valve actuators and their corresponding degrees of rotation are listed in [Table 2](#).

The setpoint for the MFC is shown in [Appendix II](#).

# 7.0 DAILY OPERATION

## 7.1 Initial Setup

- 7.1.1 Canister samples are connected to the autosampler using appropriate tubing and fittings. The sample canister valves are opened and the canister pressure gauge is monitored to assure a leak-free connection.
- 7.1.2 The GC system configuration profile, sample file list (.smp), sequence list (.seq) and method file (.mth) are set up on the Star Workstation. [Appendix I](#) has further details.
- 7.1.3 The data acquisition sequence list (.seq) and method file (.mth) are set up on the PE Data Station. [Appendix III](#) has further details.
- 7.1.4 The sample flow rate setting is confirmed on the MFC's Control/Digital Readout module. The sample volume is determined as the product of the trapping time, in minutes, times the flow rate, in  $\text{cm}^3/\text{min}$ , set on the MFC. For example:

Trapping Time:	8.0 minutes
Flow Rate:	$50.0 \text{ cm}^3/\text{min}$
Volume:	$8.00 \text{ min} \times 50.0 \text{ cm}^3/\text{min} = 400 \text{ cm}^3$

Confirmation of the actual flow rate can be done with an external flow meter.

## 7.2 Sample Preconcentration and Analysis

- 7.2.1 Samples are introduced onto the cryotrap under control of the Pre-Injection Program (PIP), an automated relay time program setup in the sample list on the Star Chromatography Workstation. The PIP valve configurations and events are shown in [Figure 1](#) through [Figure 5](#). The PIP times, relay # and status, and events are set in the sample list file (.smp) as shown in [Table 3: Pre-Injection Program Times, Relay #s, and Status](#).
- 7.2.2 At the end of the PIP, control of the analysis is passed to the GC Program. It controls heating of the cryotrap through the GC Injector Program, and the subsequent direct transfer of the trapped sample onto the GC column through the GC Relay Program. The GC Relay Programs valve configurations and events are shown in [Figure 6](#) through [Figure 8](#). The GC Relay Program times, relay # and status, and events are set in the method file (.mth) as shown in [Table 4: GC Relay Program Times, Relay #s, and Status](#).

## 7.3 Samples

- 7.3.1 A system blank is analyzed prior to calibration standards, controls and samples to evaluate the system cleanliness. System blanks are also be run after samples which contains high concentrations (>4 ppb of 1,3-butadiene or >20 ppb of benzene) to detect and eliminate possible carry-over.
- 7.3.2 Daily calibration standards are analyzed after the system blank, prior to controls or samples.
- 7.3.3 A control sample is analyzed after the system blank and calibration standards, prior to ambient air samples, in order to evaluate the accuracy of the calibration and the overall performance of the system, as shown in Table 4. The control sample contains only 1,3-butadiene, vinyl chloride, and benzene in N<sub>2</sub>. The daily results of these control standards must fall within the control limits as defined in the Method Quality Control Report.
- 7.3.4 Ambient samples are analyzed using the same sample volume as used for the calibration standard and control sample. A smaller volume is analyzed for samples containing concentrations of 1,3-butadiene or benzene that exceed the calibrated range of the analy-

sis. Smaller volumes are obtained by reducing the pre-concentration time while keeping the MFC setpoint constant.

- 7.3.5 Duplicate analyses are performed on 10% of all ambient samples analyzed. The relative percent difference of the duplicate analyses, for samples with 1,3-butadiene and/or benzene concentrations greater than 5 X Limit of Detection (LOD) are recorded and included in the method quality control report.

## 8.0 DATA ANALYSIS

- 8.1 After data acquisition, the raw data files (.raw) collected on the PE Data Station are processed by the Turbochrom™ software to produce result files (.rst). The result files contain the integrated chromatographic peak areas and retention times.
- 8.2 Chromatographic peaks in the result files for calibration standards are qualitatively identified based on retention times and analyst experience. After analyte identification, the calibration standard result files are used to calibrate the Turbochrom™ method for both retention time and concentration. The latter is based on the peak areas and the known analyte concentration in the standards.
- 8.3 After calibration of the method, chromatographic peaks in blank, control, and ambient sample result files are qualitatively identified based on relative retention times and are quantified using the calibration curve in the method.
- 8.4 A typical chromatogram of an ambient air sample is shown in [Figure 10](#).

## 9.0 QUALITY CONTROL

### 9.1 System Blank

A system blank must be analyzed before any standard or sample is run. The system blank must not contain 1,3-butadiene nor benzene at greater than two times the LOD in order to validate any subsequently analyzed samples. Trip blanks, if available, are analyzed like samples and their results are documented and evaluated.

### 9.2 Multipoint Analysis Verification

A multipoint analysis must be performed every year to verify the precision and the calibration working range. This is done by analyzing at least five simulated concentration levels of the NIST standard in triplicate.

A multipoint calibration is also required under the following conditions:

- a.) When the column is changed.
- b.) When major maintenance is performed.
- c.) When there is a change in the matrix or a reagent.

Typical calibration working range determinations are presented in [Figure 11](#) and [Figure 12](#).

### 9.3 Limit of Detection (LOD) Verification

The LOD must also be verified on a yearly basis, or when the same conditions as listed under the multipoint calibration verification occur. This is done by analyzing a minimum of seven (7) replicates of the lowest standard concentration used for calibration. This concentration must be no more than five (5) times the estimated LOD.

The LOD is estimated using the following equation, as specified in Section 9.1 of the ELB Quality Control Manual

$$\text{LOD} = t_{(n-1, 1-\alpha = 99\%)} \times S$$

$$t_{(n-1, 1-\alpha = 99\%)} = \text{student's T-distribution value at } n-1 \text{ degrees of freedom}$$

$$S = \text{standard deviation}$$

The LOD determinations are presented in [Table 5](#).

The published LODs for this method are:

1,3-butadiene .....	0.04 ppbv
benzene .....	0.2 ppbv

All subsequent LOD verifications must be equal to or less than these values.

### 9.4 Daily Calibration

A multi-point 2<sup>nd</sup> order calibration curve, forced through zero, is calculated for each analytical run from the calibration standard analyses performed during that run. A minimum of three standard concentrations are necessary to construct this type of standard curve. All calculations to produce this curve are done by the PE Turbochrom™ Data Analysis software. The Turbochrom™

method is updated after every run with the new calibration information. The method, calibration curves, the equation coefficients, and the correlation coefficients can be printed for a hardcopy record.

To be acceptable, the 2<sup>nd</sup> order correlation coefficient must be 0.98 or better.

## 9.5 Control Standard

A standard mix, containing 1,3-butadiene and benzene at concentrations within the calibration range of the GC System, is chosen to be the control standard. Analysis results of 1,3-butadiene and benzene in this standard are recorded and used to generate control charts. At least 20 data points are needed for the initial set of control limits, and any subsequent adjustment of these limits. Typical Control Charts are shown in [Figure 13](#).

The upper and lower control limits are set at  $\pm$  three times the standard deviations (SD) from the average. The upper and lower warning limits are set at  $\pm$  two times the SD from the average. If the % relative standard deviation (%RSD) of the control set is less than 5%, the %RSD is set to 5% and the SD is revised accordingly. The control are then set using the revised SD. This provides minimum upper and lower control limits of  $\pm$  15 % and upper and lower warning limits of  $\pm$  10 %. A typical dataset used for calculating control limits is given in [Table 6](#).

Control standard results must be within the established control limits for sample analyses to be valid.

## 9.6 Method Precision

Sample precision is measured by the analysis of ambient duplicate samples and the analysis of ambient collocated samples. The frequency of duplicate analyses is 10% of the total ambient samples analyzed. Control limits for the percent difference (PD) of the duplicate sample analyses is the same as the control limits for the Control Standard.

The PD for collocated sample analyses is used to evaluate method precision including both the sampling and analysis. The PD for co-located sample analyses should be within  $\pm$  25%.



## 9.7 Method Accuracy

Method accuracy is assessed by periodically analyzing other standard reference materials, if available. The results of replicate analysis of these materials must be within  $\pm 20\%$  of their reported values.

Audit samples, provided by either MLD's QM&OSB's Quality Assurance Section, or another outside source are analyzed when available. They can also be used for assessing method accuracy.

**Table 1: Function of Valves for Varian GC 3400CX System**

<b>Valve #</b>	<b>Function</b>	<b>Relay Event</b>	<b>Rotation</b>
1	Int. Std/Spike Flow Off	- 1	On
	Int. Std/Spike Flow On	+1	Off
2	Trap In Series	-4	Counterclockwise
	Trap Isolated	+4	Clockwise
3	Purge Lines – Sample Flow Off Inject Internal Standard/Spike	- 3	Counterclockwise
	Sample Flow On Load Internal Standard/Spike	+3	Clockwise
4	Inject Trap	- 2	Counterclockwise
	Load Trap	+2	Clockwise

**Table 2. List of Valve Actuators for Varian GC 3400CX System**

<b>Relay #</b>	<b># of Valve Ports</b>	<b>Rotation</b>	<b>Description</b>
2	4	90°	Series/Bypass
3	10	36°	GSV Sample on/off
4	6	60°	GSV

GSV: Gas Sampling Valve; Sx: Sample

**Table 3: Pre-Injection Program Times, Relay #s, and Status**

Time (minutes)	Relay # & Status	Events
0.00	-1-2-3-4	All valves are off (-). The sample flow is blocked and N <sub>2</sub> purge gas is allowed to flow through the MFC to the vent. The carrier gas is directed through the cryotrap to the column.
0.01	-1-2+3-4	Valve #3 is turned on (+3). This allows the sample to flow through the mass flow controller to the vent and purge the lines with new sample. The carrier gas is directed through the cryotrap to the column.
1.00	-1+2+3-4	Valve #4 is turned on (+2) and Valve #3 remains on (+3). The sample flow is directed through the cryotrap, then the MFC to vent. This starts sample loading. The carrier gas continues through the column.
9.00	-1+2-3-4	Valve #4 remains on (+2) and Valve 3 is turned off (-3). This terminates sample loading. The sample flow is blocked, allowing N <sub>2</sub> to purge the interconnecting lines and their contents onto the cryotrap. The carrier gas continues through the column.  <b><i>Note: The sample volume is varied by controlling the actions of Valve #3.</i></b>
9.10	-1+2-3+4	Valve #4 remains on (+2) and Valve #2 is turned on (+4). The cryotrap is isolated. N <sub>2</sub> purge gas is allowed to flow through the MFC to the vent. The carrier gas continues through the column.  <b><i>Note: This ends the pre-injection program events.</i></b>  <b><i><u>The activation of Valve #4 starts PE Nelson data acquisition.</u></i></b>

**Table 4: GC Relay Program Times, Relay #s, and Status**

Time (minutes)	Relay # & Status	Events
0.00	-1+2-3+4	<p>Valves #4 and #2 remain on (+2+4). The cryotrap is isolated. N<sub>2</sub> purge gas is allowed to flow through the MFC to the vent. He carrier gas continues through the column.</p> <p><b>Note: This is exactly the same state as the last step of the PIP. In fact, the last step of the PIP overwrites the first step of the GC Relay Program.</b></p>
0.02	-1-2-3+4	<p>Valve #2 remains on (+4) and Valve #4 is turned off (-2). The cryotrap is isolated and He carrier gas is directed through the column. Cryotrap pre-heating in started.</p>
1.50	-1-2-3-4	<p>Valve #2 is turned off (-4). He carrier flow is directed through the cryotrap, then to the MFC and the vent. This starts sample injection onto the column.</p> <p><b>Note: The system remains in this state throughout the run. This relay state is exactly the same as the first step in the PIP.</b></p>

**Table 5: Estimated Limits of Detection**

August 15, 1998

RUN	STANDARD	BUTADIENE			BENZENE	
		Area Counts	ppb		Area Counts	ppb
7	AAL67810	561	0.04		5719	0.25
8	AAL67810	497	0.04		5186	0.23
9	AAL67810	430	0.03		5521	0.25
10	AAL67810	481	0.04		5438	0.24
11	AAL67810	420	0.03		5222	0.23
12	AAL67810	411	0.03		5220	0.23
13	AAL67810	443	0.03		5173	0.23
	Avg.	463.29	0.03		5354.14	0.24
	Std dev	53.43	0.01		209.97	0.01
	% Rsd	11.53	15.59		3.92	4.01

3.14 x Std Dev		0.02			0.03
n =		7			7
Established LOD		0.04			0.20

**Table 6: Precision Measurements of 1,3-Butadiene and Benzene**

**Control Values for NLB57**

**Established on July 14, 1999**

**Control Gas CC113563/Calibration ALM29258**

Analysis Date		Butadiene	Benzene
July 14		0.988	3.462
		0.989	3.505
		0.982	3.538
		0.981	3.547
		0.989	3.565
		0.963	3.562
		0.980	3.530
		0.973	3.558
		0.975	3.540
		0.980	3.538
		0.963	3.531
		0.965	3.596
		0.965	3.589
		0.956	3.585
		0.957	3.605
		0.945	3.547
July 20		0.970	3.342
		0.964	3.350
		0.971	3.377
		0.968	3.551
		0.966	3.534
		0.969	3.517
July 21		0.972	3.528
		0.958	3.503
		0.961	3.511
July 26		0.973	3.288
		0.980	3.394
		0.973	3.409
		0.955	3.590
		0.956	3.569
July 28		0.980	3.462
		0.991	3.458
		0.965	3.578
Mean		0.97	3.51
Std Dev.		0.011	0.081
RSD%		1.2	2.3
PRECISION	Assigned RSD	5.0	5.0
	Assigned Std	0.049	0.175
	# of Observations	33	33
Control Limits			
	Upper Limit	1.12	4.03
	Upper Warning	1.07	3.86
	Lower Warning	0.87	3.16
	Lower Limit	0.82	2.98

**Table 7: MLD057 Standards and Controls**

<b>Date Range</b>	<b>Standard Cylinder</b>	<b>Control Cylinder</b>
01/01/96 - 09/10/96	NIST15293	C72083
09/12/96 - 11/07/96	NIST9007	C72083
11/18/96 - 03/31/97	NIST15293	C72083
04/01/97 - 09/30/97	NIST15293	AAL3113
10/01/97 - 01/12/98	AAL3113	NIST15293
01/13/98 - 02/09/98	ALM 9015	NIST15293
02/10/98 - 06/30/98	ALM 9015	ALM 9005
07/03/98 - 06/29/99	AAL3113	CC113563
07/07/99 – present	AML29258	CC113563

**Figure 1: Pre-Injection Program - Idle State**


**Time = 0.00 min - 9.10 min**

This is the standby/ready position of the GC system.

**Note:** Valves 2,3, and 4 work in combination.

### VALVE STATUS


- 1 ---- No sample is spiked
- 2 ---- Cryotrap is isolated
- 3 ---- Sample flow dead ends and N<sub>2</sub> flows through the MFC
- 4 ---- He carrier gas flows through the column





Time = 0.01 min - 9.09 min

- 1 ---- No sample is spiked
- 2 ---- Cryotrap is isolated
- +3---- Sample flows through to the MFC, purging the sample line
- 4 ---- He carrier gas flows through the column




**Figure 3: Pre-Injection Program - Sample Loaded**

Time = 1.00 min -8.10 min

### VALVE STATUS

- 1 ---- No sample spiked.
- +2---- Sample flow is directed through the cryotrap.
- +3---- Sample flows through to the MFC and the N<sub>2</sub> flow dead ends.
- 4 ---- He carrier gas flows through the column.




**Figure 4: Pre-Injection Program - Sample Loading and Cryotrap Purge**

**Time = 9.00 min -0.10 min**

# **VALVE STATUS**

- 1 ---- No sample spiked**
- +2---- Sample loading is terminated and N<sub>2</sub> flows through cryotrap**
- 3 ---- Sample flow dead ends and N<sub>2</sub> flows through to the MFC**
- 4 ---- He carrier gas flows through column**




**Figure 5: Pre-Injection Program - Injector Isolated for Heating**

Time = 9.10 min    -0.00 min

# **VALVE STATUS**

- 1 ---- No sample spiked
- +2---- Sample trap is isolated
- 3 ---- Sample flow dead ends and N<sub>2</sub> flows through the MFC
- +4---- He carrier gas flows through the column




**Figure 6: Post-Injection Program - Varian GC Control Started**

**Time = 0.00 min**

**Valve Status**

- 1 ---- No sample spiked
- +2---- Sample trap is isolated
- 3 ---- Sample flow dead ends and N<sub>2</sub> flows through the MFC
- +4---- He carrier gas flows through the column




**Figure 7: Post-Injection Program - Transition and Cryotrap Pre-Heated**

**Time = 0.02 min**

**Valve Status**

- 1 ---- No sample is spiked
- +2---- The cryotrap is isolated
- 3 ---- Sample flow dead ends and N<sub>2</sub> flows through the MFC
- 4 ---- He carrier gas flows through the column




**Figure 8: Post-Injection Program - Sample Injected onto the Column**

**Time = 1.50 min**

**Valve Status**

- 1 ---- No sample spiked
- 2 ---- He flows through the cryotrap and onto the column
- 3 ---- N<sub>2</sub> flows to the MFC
- 4 ---- He flows through the cryotrap and onto the column



**Figure 9: Injector Temperature Conversion**

Source: **ULTRA-TRACE HYDROCARBON SYSTEM**  
Operator's Manual  
Randall Bramston-Cook  
Lotus Consulting








Figure 10: Chromatograms

1. NIST 8399 Traceable Gas Mix Standard



2. Ambient Air Sample Analysis for 1,3 – Butadiene and Benzene



**Figure 11: Multipoint Analysis for 1,3-Butadiene**


July 8, 1999  
1.1 ppbv Cylinder ALM29258

	Levels of Concentration (ppb)							
Cc	25	50	75	100	200	300	400	500
PPB	0.06880	0.13750	0.20630	0.27500	0.55000	0.83000	1.10000	1.38000
1 <sup>st</sup>	3070	7130	11206	17660	39689	64649	92068	118302
2 <sup>nd</sup>	2928	7097	12243	17397	39510	62670	89940	114839
3 <sup>rd</sup>	2779	7050	11308	16530	37846	60790	88001	105331
4 <sup>th</sup>		7241	11915	16635	37260		85110	101338
Mean	2926	7092	11586	17196	39015	62703	90003	112824
Std Dev	146	40	572	591	1016	1930	2034	6716
% RSD	0.050	0.006	0.049	0.034	0.026	0.031	0.023	0.060
n	3	4	4	4	4	3	4	4

Quadratic Coefficients:  $Y=A+BX+CX^2$

Where: A=-3642.5 B=75769.6 C=5431.0

Correlation Coeff. R = 0.9982



**Figure 12: Multipoint Analysis for Benzene**

July 8, 1999

5.2 ppbv Cylinder AAL3113

	Levels of Concentration(ppb)								
cc	25	50	75	100	200	300	400	500	600
PPB	0.3250	0.6500	0.9750	1.3000	2.6000	3.9000	5.2000	6.5	7.8
1 <sup>st</sup>	19536	57384	100669	151082	351495	569932	806363	1057409	1309169
2 <sup>nd</sup>	17248	61361	103853	149208	339184	548380	789483	1018191	1262399
3 <sup>rd</sup>	18305	58042	104364	144379	330426	533419	763114	929074	1164641
4 <sup>th</sup>		59107	101298	143934	330366		732787	925484	1149486
Mean	18363	58929	102962	148223	340368	550577	786320	1001558	1245403
Std Dev	1145	2132	2002	3458	10584	18355	21797	65764	73748
% RSD	0.062	0.036	0.019	0.023	0.031	0.033	0.028	0.066	0.059
n	3	3	3	3	3	3	3	3	3

Quadratic Coefficients:  $Y = A + BX + CX^2$

Where: A=-32799 B=137341 C=3016

Correlation Coeff. R = 0.9971





Figure 13: Control Charts



## Appendix I

### Varian Star Chromatography Workstation

A Varian GC Star Workstation includes an Intel compatible PC, up to four (4) Varian Analog-to-Digital Converter Boards (ADCB), and Varian Star Chromatography software, Version 4.51. The chromatography software operates under Microsoft Windows 9.X or Microsoft Windows NT 4.0. The Star Workstation controls the gas chromatograph, documents the results of chromatographic runs, and is also capable of automating raw data collection and analysis. It is used by the ELB/OLS for automation of the Lotus Consulting/Varian instruments used in the lab. For a more detailed discussion of the Star Workstation software, including setting up methods, sequences, and sample lists, refer to the Star Chromatography Workstation Operational Manual, Varian Chromatography Systems (1993), the Varian GC Star Workstation Manual, and the Ultra-Trace Hydrocarbon System Operator's Manual, by Randall Bramston-Cook.

Although the Star Workstation is capable of data analysis and reporting, Perkin-Elmer (PE) Data Stations running PE's Turbochrom™ software is used for this purpose in the OLS. The reasons include familiarity, smaller data files, and ease in transferring to the Laboratory Information Management System (LIMS). The Star data files are also collected, however, for instrument automation/condition monitoring and as a back up for Turbochrom™ data files.

Each Gas Chromatograph (GC) serviced by the Star Workstation is assigned a separate address, as is each Analog-to-Digital Converter Board (ADCB) installed in the Workstation. Each Workstation can have up to four ADCBs and each ADCB is capable of controlling one instrument and handling two data channels. The address for each ADCB is physically set through DIP switches on the board. Each GC is assigned an address one number greater than it's associated ADCB.

For example, a Varian 3400 GC named BUTADIENE-16/17 is used to analyze butadiene and benzene. The ADCB is configured for address 16 and the attached GC is assigned address 17. It is autocontrolled by method MLD57.MTH and is attached to instrument 17. Other methods include IDLE16.MTH for system standby, BAKE16.MTH for conditioning/bakeout of the system, and PURGE16.MTH for autosampler line purging. They are used in automated sequences along with method MLD57.

A copy of the current Star Workstation method, MLD57, is included below. Although there are sections for data handling and reporting, they are not used in this analysis. They are shown in lighter type. Methods IDLE16.MTH, a stand-by method, BAKE16.MTH, a cleaning/conditioning method, are also listed.

## Varian Star Workstation Method MLD57.MTH

\*\*\*\*\*  
Star Chromatography Workstation - Method Listing  
Thu Mar 02 17:48:56 2000

Method: C:\STAR\MLD57.MTH

\*\*\*\*\*  
\*\*\*\*\*

ADC Board

\*\*\*\*\*

Module Address: 16

End Time : 24.00 minutes  
Autozero at Start : Yes  
Channel A Name : A  
Channel B Name : B  
Channel A Full Scale : 10  
Channel B Full Scale : 10

### Detector Information

-----

Detector Bunch Rate : 16 points (2.5 Hz)  
Monitor Length : 64 bunched points (25.6 seconds)  
Data File Name : mar01  
Number Files From : 1

### Integration Parameters Address 16 Channel A

-----

Subtract Blank Baseline : No  
Initial S/N Ratio : 10  
Initial Peak Width : 4 sec  
Initial Tangent Height % : 10%  
Monitor Noise : Before every run  
Measurement Type : Peak Area  
Initial Peak Reject Value : 500 counts  
Report Unidentified Peaks : No  
Report Missing Peaks : No

### Calibration Setup Address 16 Channel A

-----

Calculation Type : % (No Calibration)  
Number of Calibration Levels : 3  
Curve Origin : Include  
Curve Fit : Linear  
Weighted Regression : None)  
Replicate Treatment : Keep Replicates Separate  
Replicate Tolerance : Always add new replicates

## Varian Star Workstation Method MLD57.MTH

Out-of-Tolerance Action : No Action  
Calibration Range Tolerance : 10.0%  
Out-of-Tolerance Action : No Action

Verification Setup Address 16 Channel A  
-----

Deviation Tolerance : 100.0%  
Out-of-Tolerance Action : No Action

Peak Table Address 16 Channel A  
-----

Reference Peaks Time Windows:Width: 0.10 min. Retention Time 2.0%  
Other Peaks Time Windows :Width: 0.10 min. Retention Time 2.0%

Peak Table Empty

Time Events Table Address 16 Channel A  
-----

Time Events Table Empty

Report Format: Module ADCB Address 16 Channel A  
-----

Title : MLD51A - Butadiene and Benzene  
Print Chromatogram : No

Print Results : No  
Convert Results to ASCII? : Off

Calibration Block Reports  
Print Report : No  
Convert Report to ASCII? : Off

Print Copies : 1

\*\*\*\*\*

3400 GC

\*\*\*\*\*

Module Address: 17

GC Injector  
-----

Injector Type : Temperature Programmable  
GC Injector Oven On? : Yes  
Initial GC Injector Temperature : -50 degrees C  
Initial GC Injector Hold Time : 0.00 minutes

## Varian Star Workstation Method MLD57.MTH

Coolant To Injector/Aux. Valve On?: Yes  
Coolant Timeout : 30.00 minutes

### GC Injector Temperature Program 1

-----  
Final Temperature : 180 degrees C  
Rate : 300.0 degrees C / minute  
Hold Time : 2.00 minutes

### GC Injector Temperature Program 2

-----  
Final Temperature : 240 degrees C  
Rate : 300.0 degrees C / minute  
Hold Time : 25.04 minutes

### GC Auxiliary

-----  
Injector Type : Isothermal  
GC Auxiliary Oven On? : Yes  
GC Auxiliary Description : Valves  
Initial GC Auxiliary Temperature : 135 degrees C  
Initial GC Auxiliary Hold Time : 0.00 minutes

### GC Column

-----  
Column Oven On? : Yes  
Initial Column Temperature : 10 degrees C  
Initial Column Hold Time : 10.00 minutes  
Thermal Stabilization Time : 2.00 minutes

Coolant To Column Valve On? : Yes  
Coolant Timeout : 30.00 minutes

### GC Column Temperature Program 1

-----  
Final Temperature : 40 degrees C  
Rate : 50.0 degrees C / minute  
Hold Time : 1.00 minutes

### GC Column Temperature Program 2

-----  
Final Temperature : 120 degrees C  
Rate : 8.0 degrees C / minute  
Hold Time : 1.00 minutes

### GC Column Temperature Program 3



## Varian Star Workstation Method MLD57.MTH

-----  
Final Temperature : 220 degrees C  
Rate : 50.0 degrees C / minute  
Hold Time : 3.40minutes

### GC Column A Parameters

Installed? : Yes  
Length : 75.0 meters  
Diameter : 450 microns  
Carrier Gas : Helium

### GC Column B Parameters

Installed? : Yes  
Length : 75.0 meters  
Diameter : 450 microns  
Carrier Gas : Helium

### GC Detector Heater

-----  
Detector Heater On? : Yes  
Detector Temperature : 245 degrees C

### GC Detector A

-----  
Detector Type : PID  
Detector On? : Yes  
Attenuation : 8  
Detector Range : 12  
Autozero at GC Ready? : No

### GC Detector B

-----  
Detector Type : ECD  
Detector On? : No  
Attenuation : 8  
Detector Range : 10  
Autozero at GC Ready? : Yes

### Autosampler

-----  
Autosampler Type : 8134 SSV

### GC Relays

-----

## Varian Star Workstation Method MLD57.MTH

Relay Time Program : Use  
Initial Relay States : -2-4  
Relay Initial Conditions at Run End?: No

Relay Program 1

-----

Relay Time: 0.02 State ---4

Relay Program 2

-----

Relay Time : 1.50 State ----

GC Stripchart

-----

Stripchart On? : No

## Varian Star Workstation Method IDLE16.MTH

\*\*\*\*\*  
Star Chromatography Workstation  
Method Listing Fri Mar 03 17:03:53 2000

Method: C:\STAR\IDLE16.MTH

\*\*\*\*\*  
\*\*\*\*\*

ADC Board

\*\*\*\*\*

Module Address: 16

End Time : 30.00 minutes  
Autozero at Start : Yes  
Channel A Name : PID  
Channel B Name : ECD  
Channel A Full Scale : 10  
Channel B Full Scale : 10

### Detector Information

-----

Detector Bunch Rate : 2 points (20.0 Hz)  
Monitor Length : 64 bunched points (3.2 seconds)  
Data File Name : idle2  
Number Files From : 1

### Integration Parameters Address 16 Channel A

-----

Subtract Blank Baseline : No  
Initial S/N Ratio : 5  
Initial Peak Width : 4 sec  
Initial Tangent Height % : 10%  
Monitor Noise : Before every run  
Measurement Type : Peak Area  
Initial Peak Reject Value : 1000 counts  
Report Unidentified Peaks : Yes  
Report Missing Peaks : No

### Calibration Setup Address 16 Channel A

-----

Calculation Type : % (No Calibration)  
Number of Calibration Levels : 10  
Curve Origin : Force  
Curve Fit : Linear  
Weighted Regression : (None)  
Replicate Treatment : Keep Replicates Separate  
Replicate Tolerance : Always add new replicates

## Varian Star Workstation Method IDLE16.MTH

Out-of-Tolerance Action : No Action  
Calibration Range Tolerance : 10.0%  
Out-of-Tolerance Action : No Action

Verification Setup Address 16 Channel A

-----  
Deviation Tolerance : 100.0%  
Out-of-Tolerance Action : No Action

Peak Table Address 16 Channel A

-----  
Reference Peaks Time Windows:Width: 0.10 min. Retention Time 2.0%  
Other Peaks Time Windows :Width: 0.10 min. Retention Time 2.0%

Peak Table Empty

Time Events Table Address 16 Channel A

-----  
Time Events Table Empty

Report Format: Module ADCB Address 16 Channel A

-----  
Title :  
Print Chromatogram : Yes  
Chromatogram Options :  
Start Retention Time : 0.00 minutes  
End Retention Time: 1440.00 minutes  
Length in Pages : 1  
Initial Chart Speed : 0.0 cm/min  
Minutes per Tick : 1.0  
Autoscale : On  
Time Events : On  
Chromatogram Events : On  
Retention Times : On  
Peak Names : On  
Baseline : On  
Print Results : Yes  
Results Options:  
Units :  
Number of Decimal Digits : 4  
Show Peak Group Totals : No  
Run Log : Off  
Error Log : Off  
Calibration Report : Off  
Revision Log : On  
Notes : Off

## Varian Star Workstation Method IDLE16.MTH

Convert Results to ASCII? : Off

### Calibration Block Reports

Print Report : No

Convert Report to ASCII? : Off

Print Copies : 1

\*\*\*\*\*

3400 GC

\*\*\*\*\*

Module Address: 17

GC Injector

Injector Type : Temperature Programmable

GC Injector Oven On? : Yes

Initial GC Injector Temperature : 200 degrees C04444

Initial GC Injector Hold Time : INFINITE

Coolant To Injector/Aux Valve On?: No

Coolant Timeout : INFINITE

GC Auxiliary

-----

Injector Type : Isothermal

GC Auxiliary Oven On? : Yes

GC Auxiliary Description : Valving

Initial GC Auxiliary Temperature : 135 degrees C

Initial GC Auxiliary Hold Time : 0.00 minutes

GC Column

-----

Column Oven On? : Yes

Initial Column Temperature : 100 degrees C

Initial Column Hold Time : INFINITE

Thermal Stabilization Time : 2.00 minutes

Coolant To Column Valve On? : No

Coolant Timeout : INFINITE

GC Column A Parameters

-----

Installed? : Yes

Length : 75.0 meters

Diameter : 450 microns

Carrier Gas : Helium

## Varian Star Workstation Method IDLE16.MTH

### GC Column B Parameters

-----

Installed? : No  
Length : 0.0 meters  
Diameter : 0 microns  
Carrier Gas : Helium

### GC Detector Heater

-----

Detector Heater On? : Yes  
Detector Temperature : 245 degrees C

### GC Detector A

-----

Detector Type : PID  
Detector On? : Yes  
Attenuation : 8  
Detector Range : 10  
Autozero at GC Ready? : Yes

### GC Detector B

-----

Detector Type : Not Used  
Detector On? : No

### Autosampler

-----

Autosampler Type : 8134 SSV

### GC Relays

-----

Relay Time Program : Use  
Initial Relay States : -2--  
Relay Initial Conditions at Run End? : Yes

### GC Stripchart

-----

Stripchart On? : No

## Varian Star Workstation Method BAKEOU16.MTH

\*\*\*\*\*

Star Chromatography Workstation

Method Listing Fri Mar 03 17:04:57 2000

Method: C:\STAR\BAKEOU16.MTH

\*\*\*\*\*

\*\*\*\*\*

ADC Board

\*\*\*\*\*

Module Address: 16

End Time : 20.20 minutes

Autozero at Start : Yes

Channel A Name : A

Channel B Name : B

Channel A Full Scale :

Channel B Full Scale :

Detector Information

-----

Detector Bunch Rate : 16 points (2.5 Hz)

Monitor Length : 64 bunched points (25.6 seconds)

Data File Name : star

Number Files From : 1

Integration Parameters Address 16 Channel A

-----

Subtract Blank Baseline : No

Initial S/N Ratio : 5

Initial Peak Width : 4 sec

Initial Tangent Height % : 10%

Monitor Noise : Before every run

Measurement Type : Peak Area

Initial Peak Reject Value : 1000 counts

Report Unidentified Peaks : Yes

Report Missing Peaks : No

Calibration Setup Address 16 Channel A

-----

Calculation Type : % (No Calibration)

Number of Calibration Levels : 9

Curve Origin : Force

Curve Fit : Linear

Weighted Regression : (None)

Replicate Treatment : Keep Replicates Separate

Replicate Tolerance : Always add new replicates

## Varian Star Workstation Method BAKEOU16.MTH

Out-of-Tolerance Action : No Action  
Calibration Range Tolerance : 10.0%  
Out-of-Tolerance Action : No Action

Verification Setup Address 16 Channel A

-----  
Deviation Tolerance : 100.0%  
Out-of-Tolerance Action : No Action

Peak Table Address 16 Channel A

-----  
Reference Peaks Time Windows:Width: 0.10 min. Retention Time 2.0%  
Other Peaks Time Windows :Width: 0.10 min. Retention Time 2.0%

Peak Table Empty

Time Events Table Address 16 Channel A

-----  
Time Events Table Empty

Report Format: Module ADCB Address 16 Channel A

-----  
Title : Bakeout Chromatogram - Butadiene  
Method 57

Print Chromatogram : Yes  
Chromatogram Options:  
Start Retention Time : 0.00 minutes  
End Retention Time : 1440.00 minutes  
Length in Pages : 1  
Initial Chart Speed : 0.0 cm/min  
Minutes per Tick : 1.0  
Autoscale : On  
Time Events : On  
Chromatogram Events : On  
Retention Times : On  
Peak Names : On  
Baseline : On

Print Results : No  
Convert Results to ASCII? : Off

Calibration Block Reports

Print Report : No  
Convert Report to ASCII? : Off



## Varian Star Workstation Method BAKEOU16.MTH

Print Copies : 1

\*\*\*\*\*

3400 GC

\*\*\*\*\*

Module Address: 17

GC Injector

-----

Injector Type : Temperature Programmable

GC Injector Oven On? : Yes

Initial GC Injector Temperature : 200 degrees C

Initial GC Injector Hold Time : 0.50 minutes

Coolant To Injector/Aux Valve On? : Yes

Coolant Timeout : INFINITE

GC Injector Temperature Program 1

-----

Final Temperature : 260 degrees C

Rate : 300.0 degrees C / minute

Hold Time : 19.50 minutes

GC Auxiliary

-----

Injector Type : Isothermal

GC Auxiliary Oven On? : Yes

GC Auxiliary Description : Valving

Initial GC Auxiliary Temperature : 135 degrees C

Initial GC Auxiliary Hold Time : 0.00 minutes

GC Column

-----

Column Oven On? : Yes

Initial Column Temperature : 190 degrees C

Initial Column Hold Time : 0.00 minutes

Thermal Stabilization Time : 2.00 minutes

Coolant To Column Valve On? : Yes

Coolant Timeout : INFINITE

## Varian Star Workstation Method BAKEOU16.MTH

### GC Column Temperature Program 1

-----  
Final Temperature : 200 degrees C  
Rate : 50.0 degrees C / minute  
Hold Time : 20.00 minutes

### GC Column A Parameters

-----  
Installed? : No  
Length : 0.0 meters  
Diameter : 0 microns  
Carrier Gas : Helium

### GC Column B Parameters

-----  
Installed? : No  
Length : 0.0 meters  
Diameter : 0 microns  
Carrier Gas : Helium

### GC Detector Heater

-----  
Detector Heater On? : Yes  
Detector Temperature : 245 degrees C

### GC Detector A

-----  
Detector Type : PID  
Detector On? : Yes  
Attenuation : 8  
Detector Range : 10  
Autozero at GC Ready? : No

### GC Detector B

-----  
Detector Type : ECD  
Detector On? : No  
Attenuation : 8  
Detector Range : 1  
Autozero at GC Ready? : No

### Autosampler

-----  
Autosampler Type : 8134 SSV

## Varian Star Workstation Method BAKEOU16.MTH

GC Relays

-----

Relay Time Program : Use

Initial Relay States : ----

Relay Initial Conditions at Run End? : No

GC Stripchart

-----

Stripchart On? : No

---

## APPENDIX II

### Additional Setpoints

GC Column Stand by Temperature-----150 °C

#### He Carrier Gas:

Set to 5.00 cm<sup>3</sup>/minute

Flow Controller-----Set digital gauge at 2.00 (~ cm<sup>3</sup>/minute)

Pressure Regulator-----Set analog gauge to 49.8 psi

#### N<sub>2</sub> Make-up Gas:

PID ----- Set N<sub>2</sub> flow to 13.00 cm<sup>3</sup>/minute

ECD -----Set N<sub>2</sub> flow to 17.00 cm<sup>3</sup>/minute

The ECD is in series with the PID. It is not used in this method.

**Total He and N<sub>2</sub> Gas Flow** -----35.00 cc/minutes

#### N<sub>2</sub> Purge Gas:

Flow Controller-----Set digital gauge at 35.0 (~ cm<sup>3</sup>/minute)

N<sub>2</sub> purge Pressure Gauge -----Reads > 60 psi

#### Mass Flow Controller (MFC):

Set sampling flow rate to 50 cm<sup>3</sup>/minute

Set -----49.9% of full scale

Read -----50.0% of full scale

-----100 cm<sup>3</sup>/minute equals 100% full scale

#### Required Regulator Pressures:

He - Carrier Gas-----80 psi

N<sub>2</sub> - Valve Actuator Gas -----80 psi

N<sub>2</sub> - Make-up Gas-----80 psi

N<sub>2</sub> – Dryer Purge Gas -----80 psi

---

## APPENDIX III

### Perkin-Elmer Nelson Turbochrom Data System

#### Data Station

The Perkin-Elmer (PE) Data Station includes an Intel based PC and PE-Nelson 2700/Turbochrom™ Chromatography software, and one or more PE-Nelson Model 970 Analog/Digital Converters. Up to 15 A/D converters can be attached to a single Data Station. In practice, only two (2) to four (4) A/Ds are used per Data Station. Each A/D can support two synchronous analog data channels.

The PE-Nelson 2700/Turbochrom™ Chromatography software, Version 4.1.2 <2F12>, operates under Microsoft Windows 9.x, but will not function under the Windows NT. For a more detailed discussion of the Turbochrom™ software, please refer to the Turbochrom User's Guide, Volumes 1 and 2.

#### File Naming Conventions

Turbochrom™ uses the first four (4) characters to form a data files base name. A three-digit cycle number is then appended to the selected base name and the data files are stored in a specified directory. The same base name is used for storing raw, result, baseline, and modified raw files, but the three-character cycle number differentiates one file from another. If a duplicate file name exists, the system will rename the current sample by appending letters. This will cause problems with the eventual transfer of data to LIMS. To prevent this and to allow easy identification of data, the following naming convention is used:

1. Use only a 4 character base name for data runs.
2. The first two characters are the month code.
3. The third and fourth characters are the date code
4. The fifth through seventh characters are automatically assigned by the data system.

The applicable month codes are:

January	=	JA	May	=	MY	September	=	SE
February	=	FB	June	=	JN	October	=	OC
March	=	MR	July	=	JL	November	=	NV
April	=	AP	August	=	AG	December	=	DC1

---

For example:

<b>FB</b>	<b>01</b>	<b>020</b>	<b>=</b>	<b>FEBRUARY 1st , RUN #20</b>
Month Code	Day Code	Cycle No.		

### Sample Naming Conventions

Within a sequence or during a manual download, a sample name this must be entered. The sample name must be uniquely formatted for LIMS to recognize.

- Enter a system blank as: **Blank**
- Enter a Calibration standard as: **ALM29258 or XLM20258**
- Enter a control sample as: **C113563 or XC113563**

Ambient air samples are assigned an 8-character bar code number. When logged into the LIMS, it receives a 9-digit number that starts with a 2. For Toxics, the bar code number begins with TX, followed by 6 digits. The sample name, for Turbochrom™, includes the bar code number, plus a site code. A space is used to separate the bar code section used for LIMS and the site information. A typical sample name would be entered as:

**TX001259 MX**

A duplicate, or replicate analyses, receives an additional alpha character, A through J, at the end of the bar code number. A through J represents the 2<sup>nd</sup> through 10<sup>th</sup> replicates. The duplicate for the above sample would be entered as:

**TX001259A MX**

### Turbochrom Set-Up

The method file, MLD57A.MET follows.

---

Turbochrom Method File: C:\TC4-2F12\MLD57.MTH  
Created by : HML on : 11/6/98 10:11 AM  
Edited by : Ferry/ FNN on : 3/2/00 10:20 AM  
Description :

Number of Times Edited : 82  
Number of Times Calibrated : 2407

Instrument Conditions :  
7/3 Switched relay 2 & 4. (EMcC/MRP)  
MLD51A.MTH runtime 30 min, multiple level cals,  
new cal std, new baseline events

Instrument Control Method:  
Instrument name : 970\_BUTADIENE

Interface Parameters :  
Delay Time : 0.00 min.  
Run Time : 30.00 min.  
Sampling Rate : 2.0000 pts/s  
Interface Type : 900  
Analog Voltage Input : 1000 mV  
Data will be collected from channel A

Timed Events:  
There are no timed events in the method

Real Time Plot Parameters :  
Channel A -- Pages: 1 Offset: 0.000 mV Scale: 10.000 mV  
Channel B -- Pages: 1 Offset: 0.000 mV Scale: 1000.000 mV

Processing Parameters :  
Bunch Factor : 1 points  
Noise Threshold : 1 uV  
Area Threshold : 25.00 uV  
Peak Separation Criteria  
Width Ratio : 0.200  
Valley-to-Peak Ratio : 0.010  
Exponential Skim Criteria  
Peak Height Ratio : 5.000  
Adjusted Height Ratio : 4.000  
Valley Height Ratio : 3.000

Baseline Timed Events :  
Event #1 - -P at 0.100  
Event #2 - +P at 7.150  
Event #3 - S at 8.047  
Event #4 - +X at 8.204

---

Event #5 - S at 8.447  
Event #6 - -P at 8.500  
Event #7 - +P at 19.000  
Event #8 - S at 19.285  
Event #9 - +X at 19.553  
Event #10 - S at 19.885  
Event #11 --P at 19.945

Annotated Replot Parameters :

Offset will be autozeroed  
Plot Scale : 35.000 mV  
Number of Pages : 1  
Plot Title : Chromatogram  
X-Axis Label : Time [min]  
Y-Axis Label : Response [mV]  
Orientation : Landscape  
Retention Labels : Top of Plot  
Component Labels : Actual Time  
Start Time : 0.00  
End Time : 30.00

Report Format files :

No report format files given

User Programs :

No user programs will be executed

Global Information :

Default Sample Volume : 1.000 ml  
Quantitation Units : ppb  
Void Time : 0.000 min  
Correct amounts during calibration : YES  
Reject outliers during calibration : NO  
An External Standard calibration will be used  
Unknown peaks will use the response factor of the nearest component

Component Information :

Vinyl Chloride  
Component Type : Single Peak Component  
Retention Time : 7.571 min  
Search Window: 0.00 s, 1.30 %  
Reference Component:  
Find peak closest to expected RT in window  
Calibrating Area versus Amount using a 2nd Order Fit  
Curve will be forced through the origin  
Amounts will not be scaled prior to the regression  
Weighting factor for the regression: None



---

User Values:

Label :  
Value 1: 0.000000  
Value 2: 0.000000  
Value 3: 0.000000  
Value 4: 0.000000  
Value 5: 0.000000

Calibration Levels:

Level Name	Amount	Area	Height	ISTD Re
50	0.2050	739.50	175.04	-----
100	0.4100	1437.00	356.00	-----
400	1.6400	6750.41	1713.01	-----

Calibration Curve:  $y = (0.000000) + (3360.531646)x + (460.471100)x^2$   
R-squared : 0.999939

Butadiene

Component Type : Single Peak Component  
Retention Time : 8.177 min  
Search Window: 0.00 s, 2.00 %

Reference Component:

Find peak closest to expected RT in window

Calibrating Area versus Amount using a 2nd Order Fit

Curve will be forced through the origin

Amounts will not be scaled prior to the regression

Weighting factor for the regression: None

User Values:

Label :  
Value 1: 0.000000  
Value 2: 0.000000  
Value 3: 0.000000  
Value 4: 0.000000  
Value 5: 0.000000

Calibration Levels:

Level Name	Amount	Area	Height	ISTD Re
50	0.1380	2182.00	473.60	-----
100	0.2750	5062.50	1096.04	-----
400	1.1000	26932.00	5785.37	-----
500	1.3800	34776.00	7492.71	-----

Calibration Curve:  $y = (0.000000) + (18167.096115)x + (5275.56425)x^2$   
R-squared: 0.999064

Benzene

Component Type : Single Peak Component  
Retention Time : 19.390 min  
Search Window: 0.00 s, 1.50 %

Reference Component:  
Find largest peak in window  
Calibrating Area versus Amount using a 2nd Order Fit  
Curve will be forced through the origin  
Amounts will not be scaled prior to the regression  
Weighting factor for the regression: None  
User Values:

Label :  
Value 1: 0.000000  
Value 2: 0.000000  
Value 3: 0.000000  
Value 4: 0.000000  
Value 5: 0.000000

Calibration Levels:

Level	Name	Amount	Area	Height	ISTD Re
50		0.6500	18766.50	5379.35	-----
100		1.3000	41897.75	13079.19	-----
400		5.2000	222028.88	70972.46	-----
500		6.5000	283312.50	91407.00	-----

Calibration Curve:  $y = (0.000000) + (32437.864321)x + (1787.39247)x^2$   
R-squared : 0.998938

Calibration Replicate Lists:

Component: Vinyl Chloride

Level : 50

Area	Height	Vol	Adj	Amt	ISTD	Response	ISTD	Amount
Date/Time								

739.50	175.04	0.2050	-----	-----	3/2/00
--------	--------	--------	-------	-------	--------

Level : 100

Area	Height	Vol	Adj	Amt	ISTD	Response	ISTD	Amount
Date/Time								

1437.00	356.00	0.4100	-----	-----	3/2/00
---------	--------	--------	-------	-------	--------

Level : 400

Area	Height	Vol	Adj	Amt	ISTD	Response	ISTD	Amount
Date/Time								

6750.41	1713.01	1.6400	-----	-----	3/2/00
---------	---------	--------	-------	-------	--------

Component: Butadiene

Level : 50

Area	Height	Vol	Adj	Amt	ISTD	Response	ISTD	Amount
Date/Time								

---

2182.00	473.60	0.1380	-----	-----	3/2/00
Level : 100					
Area	Height	Vol Adj Amt	ISTD Response	ISTD Amount	
Date/Time					

---

5062.50	1096.04	0.2750	-----	-----	3/2/00
Level : 400					
Area	Height	Vol Adj Amt	ISTD Response	ISTD Amount	
Date/Time					

---

26932.00	5785.37	1.1000	-----	-----	3/2/00
Level : 500					
Area	Height	Vol Adj Amt	ISTD Response	ISTD Amount	
Date/Time					

---

34776.00	7492.71	1.3800	-----	-----	3/2/00
Component: Benzene					
Level : 50					
Area	Height	Vol Adj Amt	ISTD Response	ISTD Amount	
Date/Time					

---

18766.50	5379.35	0.6500	-----	-----	3/2/00
Level : 100					
Area	Height	Vol Adj Amt	ISTD Response	ISTD Amount	
Date/Time					

---

41897.75	13079.19	1.3000	-----	-----	3/2/00
Level : 400					
Area	Height	Vol Adj Amt	ISTD Response	ISTD Amount	
Date/Time					

---

222028.88	70972.46	5.2000	-----	-----	3/2/00
Level : 500					
Area	Height	Vol Adj Amt	ISTD Response	ISTD Amount	
Date/Time					

---

283312.50	91407.00	6.5000	-----	-----	3/2/00
-----------	----------	--------	-------	-------	--------

---