California Environmental Protection Agency

Air Resources Board

PROCEDURE FOR THE DETERMINATION OF C₂ TO C₁₂ HYDROCARBONS IN AUTOMOTIVE EXHAUST SAMPLES BY GAS CHROMATOGRAPHY

Standard Operating Procedure No. MLD 102 / 103 Revision 2.0

September 1996

Organic Analysis Section Southern Laboratory Branch Monitoring & Laboratory Division 9528 Telstar Avenue El Monte, CA 91731

This procedure has been reviewed by the staff of the Air Resources Board and approved for publication. Approval does not signify that the contents necessarily reflect the views and policies of the Air Resources Board, nor does the mention of trade names or commercial products constitute endorsement or recommendation for use.

1 Introduction

1.1 This document describes a method of analyzing, by gas chromatography (GC), C₂ to C₅ (light-end) hydrocarbons and C₆ to C₁₂ (mid-range) hydrocarbons, in the range of 20 to 10000 parts per billion carbon (ppbC), from automotive source samples. It is derived from the Air Resources Board Methods 1002 and 1003 (Ref. 9.1). This method does not include sample collection procedures.

2 Method Summary

- 2.1 This procedure uses two gas chromatographs (GCs), one (the light-end GC) for analyzing C_2 to C_5 hydrocarbons and the other (the mid-range GC) for analyzing C_6 to C_{12} hydrocarbons, each with its own automated cryogenic preconcentration system (cryotrap).
- 2.2 For routine motor vehicle exhaust testing, the vehicle is tested according to the Federal Test Procedure (FTP, Ref. 9.2), using a dynamometer (dyno) and constant volume sampler (CVS) to dilute the exhaust for sampling.
- 2.3 Samples are also received from CVS testing using non-FTP driving cycles, Sealed Housing Evaporative Determinations (SHEDs, Ref. 9.2), gas standard cylinders or canisters, and hydrocarbon-containing samples from other miscellaneous sources.
- 2.4 The samples are received by the laboratory in Tedlar bags or in stainless steel canisters. The sample is injected into the GC by means of the automatic gas phase sampling valve system. Separation of the sample hydrocarbon mixture into its components takes place in the chromatographic column. The flame ionization detector (FID) is used for detection and quantification.
- 2.5 The hydrocarbon concentrations are determined by integrating the peak areas and using response factors determined from National Institute of Standards and Technology (NIST)-traceable standards during calibration.
- 2.6 The computerized GC data acquisition system identifies the hydrocarbon associated with each chromatographic peak (Ref. 9.3)

3 Interferences and Limitations

3.1 Any component (e.g., Freon group, Tedlar bag contaminants) present in the sample having a retention time similar to that of any hydrocarbon being measured is an interferent. Therefore, proof of chemical identity requires confirmation by other methods and instrumentation, e.g., gas chromatography/mass spectrometry (GC/MS).

- 3.2 To maximize sample integrity, sample bags should not leak or be exposed to bright light or excessive heat. Sampling bags must be shielded from direct sunlight to avoid reactions due to reactive hydrocarbons. The compound 1,3-butadiene (a light-end hydrocarbon), most of which is in cold-start exhaust test bag no. 1, is unstable. Therefore, all cold-start bag no. 1 samples must be analyzed for the light-end hydrocarbons within 8 hours. Cold-start bag no. 2, bag no. 3 and background must be analyzed within 24 hours, although analysis within 8 hours is recommended.
- 3.3 The light-end analysis is calibrated with propane. The propane response factor is used for all light-end hydrocarbons. This procedure will introduce an error if the propane response is not representative of the responses of all other light-end compounds.
- 3.4 The mid-range analysis is calibrated with benzene or propane (Section 5.2.2.1). The benzene (or propane) response factor is used for all mid-range hydrocarbons and ethers, such as MTBE. This procedure will introduce an error if the benzene (or propane) response is not representative of the responses of all other mid-range compounds.

4 Instrumentation and Apparatus

- 4.1 Tedlar bags, 2 mils in thickness, nominally 5 to 10 liters in capacity and equipped with Swagelok QuikConnect fittings, are used to contain the samples from the dyno test facility.
- 4.2 Two gas chromatographs, VARIAN 3400 or 3600 or equivalent, equipped with gas sampling valve systems (as described in Ref. 9.3) and FIDs, are connected together by their sampling lines and electronic start signals.
- 4.3 The light-end GC uses an analytical column, Al203/KCl PLOT [50 meters (m) x 0.32 millimeters (mm) inner diameter (ID)] and, to prevent water damage to the PLOT column, a Carbowax WCOT (25 m x 0.53 mm ID) precolumn; the mid-range GC uses a DB-1 WCOT analytical column [40 m x 0.18 mm ID with a 0.4 micron (u) film thickness or 60 m x 0.32 mm ID with a 1 u film thickness]. Precolumns and analytical columns are purchased from CHROMPACK, Inc., J&W or equivalent.
- 4.4 A cryotrap of nickel tubing [10" x 1/16" outer diameter (OD), 0.040" ID] is used to concentrate the sample.
- 4.5 A PC-based data acquisition system (Varian Star or equivalent) is used to quantitate peak areas.

5 Reagents and Materials

- 5.1 Gases used to support the GC analysis shall be of the following purities:
- 5.1.1 Helium shall have a minimum purity of 99.9999% and pass through a molecular sieve trap to ensure purity before entering the GC system.
- 5.1.2 Hydrogen shall have a minimum purity of 99.999%.
- 5.1.3 Zero nitrogen shall have a minimum purity of 99.998%.
- 5.1.4 Ultra-zero air shall have less than 0.1 part per million (ppm) hydrocarbon (as methane).
- Gas standards with concentrations in the range of 100 to 10000 ppbC are used for daily calibration, daily quality control and linearity/limit of detection (LOD) determinations.
- 5.2.1 Primary calibration gas standards, purchased from NIST, are used to verify the concentrations of working standards used for daily calibration.
- 5.2.2 The working standard is purchased from Scott Specialty or other vendor.
 - 5.2.2.1 The working calibration standard contains at least propane (approximately 3000 ppbC) for the light-end calibration analysis and benzene (approximately 1500 ppbC) for the mid-range calibration analysis.
 - 5.2.2.2 If propane is shown to have a mid-range response factor that is representative of a broad range of mid-range compounds, propane may be used to calibrate both the light-end and mid-range analyses. In that case, the calibration standard would be required to contain only propane.
 - 5.2.2.3 The concentration(s) of the calibration compound(s) is verified (\pm 2%) by comparison with the NIST standard(s).
 - 5.2.2.3 Additional compounds may be included in the calibration standard for response factor comparisons, retention time confirmations, etc, if desired.
- 5.2.3 The quality control gas mixture is also purchased from Scott Specialty or other vendor.
 - 5.2.3.1 It contains at least the following gases: ethene, propane, n-butane, and 2-methylpropene for the light end analysis (Ref. 9.1, Method 1002, Section 5.6) and n-hexane, n-octane, n-decane, benzene, toluene, and m- or p-xylene for the mid-range analysis (Ref. 9.1, Method 1003, Section 5.6).
 - 5.2.3.2 These ten compounds are monitored for quality control purposes (Section 8).
 - 5.2.3.3 Additional compounds may be included in the control standard for response factor comparisons, retention time verifications, etc., if desired.

- 5.2.4 Table 1 shows the composition of the calibration and control standard gases in use at this date.
- 5.2.5 Figures 1 through 4 show the light-end and mid-range chromatograms of the current calibration standards.
- 5.4 Liquid nitrogen (LN2) is used to cool the cryogenic traps and column ovens.

6 **Procedure**

- 6.1 The operating conditions for the gas chromatographs are given in Table 2.
- 6.1.1 Each GC pair may have some deviation of these conditions to optimize the operation.
- 6.2 Samples in Tedlar bags (or stainless steel canisters) are attached to the gas phase autosampling valve system and simultaneously drawn through a fixed-volume sample loop in the light-end GC and through another fixed sample loop in the mid-range GC, with a vacuum pump [variable, usually at about 300 mL/min for 2 minutes (variable)]. The volume of the sample loop may vary for each GC, depending on its intended use. For the light-end GC it is usually 10 milliliters (mL), nominal, and for the mid-range GC it varies from 7 to 30 mL
- 6.3 For the light-end GC the sample is eluted through the precolumn and loaded into the cryotrap; for the mid-range GC the sample is loaded directly into the cryotrap.
- 6.4 After all the light-end hydrocarbons are eluted through the precolumn, the light-end GC precolumn is back flushed to prevent the heavy hydrocarbons and water from eluting onto the analytical (PLOT) column. Helium is flushed backward through the precolumn out to vent. The cryotraps for both GCs are isolated and heated.
- 6.5 When the cryotrap is heated for a specific time (corresponding to 200°C), the isolation valve is opened and the samples are injected into the analytical column of each GC.
- As each separated hydrocarbon elutes from the column through the FID, a signal response (peak) is generated and recorded by the data system.
- 6.7 The data system generates separate light-end and mid-range reports, each containing the chromatograms and peak identifications and corresponding concentrations.
- 6.7 The peak identifications and peak concentrations determined by the data system are reviewed and, if necessary, corrected using the following procedure and criteria:
 - (1) the relative retention indices from GC/MS analyses are used to help confirm peak

identification.

- (2) the primary peak identification is determined by the data system, using the retention times based on reference calibration runs.
- (3) the relative peak heights of the sample run ("fingerprint") are matched with the typical fingerprint seen in the past sample runs.
- (4) the relative peak retention times of the sample run are compared with those of the reference runs; and
- (5) any peak with reasonable doubt is labeled 'Unidentified'.
- 6.8 After running an extremely "dirty" sample, the analyst may run a blank or flush the sampling system with zero nitrogen before proceeding to the next sample, as there may be sample carry-over.
- 6.9 All peaks identified as target compounds at or above the LOD are reported.
- 6.9.1 Target compounds (Ref. 9.1, Appendix 1) which coelute are reported as the major component as determined by the analysis of several samples by GC/MS or other methods.
- 6.9.2 An exception to this is m- and p-xylene. It has been stipulated that two-thirds of the coeluting peak is m-xylene and one-third is due to p-xylene. This ration of 2:1 was developed from GC/MS data and fuel profiles. This ratio is not used for peak identification and reporting of NMHC; the peak is reported as "m- & p-xylene". However, this ratio is used to calculate reactivity (Ref. 9.1, Parts E and G).

7 Calculations

- 7.1 A single-point calibration is performed daily with a NIST-traceable secondary calibration standard.
- 7.2 The hydrocarbon concentrations, in ppbC, are calculated by the data system using an external standard method.

Sample Concentration = $\frac{\text{Sample Peak Area}}{\text{Response Factor}}$

| 7.5 | The response rac | to (Ri) is calculated by. |
|-----|------------------|---|
| | RF = | Peak Area of calibration standard component |
| | Con | centration of NIST-traceable calibration standard component |

7.3.1 The light-end analysis is calibrated with propane.

The response factor (RF) is calculated by:

7.3.2 The mid-range may be calibrated with either benzene or propane (Section 5.2).

8 Quality Control

73

- 8.1 A blank (pure nitrogen or helium) is run once daily before running the calibration standard, control standard, and samples. All target hydrocarbon concentrations from the blank analysis must be below the LOD before the analysis may proceed.
- 8.2 One run of the calibration standard is performed per day to generate the response factors needed for quantifying sample analyses.
- 8.3 One run of the quality control standard is performed daily.
- 8.3.1 A quality control chart is maintained for each monitored hydrocarbon compound in the control standard. The control charts, used on a daily basis, establish that the method is "incontrol". The following describes how a typical control chart is constructed:
 - (1) obtain at least 20 daily control standard results;
 - (2) calculate the control standard mean concentration and standard deviation for the monitored hydrocarbon; and
 - (3) create a control chart for the monitored hydrocarbon by placing the dates on the x-axis and the concentrations on the y-axis. Establish an upper and lower warning limit at two standard deviations (2s) above and below the average concentration. Establish an upper and lower control limit at three standard deviations (3s) above and below the average concentration.
- 8.3.2 Measurements of all ten monitored hydrocarbon compounds contained in the control standard must be within the control limits ("in-control") before sample analysis may proceed. The results of all control monitored compounds are plotted on control charts. Values which exceed three standard deviations above or below the mean are considered to be "out of control". Values which exceed two standard deviations above or below the mean on two consecutive analysis days are also considered to be out of control on the second day. If one or more monitored hydrocarbon compounds are out of control, it may be necessary to inspect and repair the GC, and rerun the calibration and/or control standards until the control standard criteria are met; when the QC criteria are met, sample analysis can continue. Figure 5 demonstrates a typical QC chart.
- 8.4 A duplicate analysis of one sample bag is performed at least once a day. The relative

percent difference (RPD) is calculated for each duplicate run:

RPD (%) = <u>Difference between duplicate and original measurements X 100</u>

Average of duplicate and original measurements

For each monitored hydrocarbon compound in the control standard, the allowable RPD depends on the average concentration of the duplicate runs, as shown in the following table (Ref. 9.1):

| Average Measurement for Duplicate Runs | Allowable RPD (%) |
|--|-------------------|
| 1 to 10 times LOD | 100 |
| 10 to 20 times LOD | 30 |
| 20 to 50 times LOD | 20 |
| Greater than 50 times LOD | 15 |

The results from duplicate analysis must meet the criteria above for all monitored hydrocarbon compounds in the control standard for Methods 1002 and 1003 (Ref.9.1). If the criteria are not met, the sample must be rerun. If the criteria are still not met, all sample results for the day from this instrument must be deleted and the samples reanalyzed. Figure 6 shows a typical duplicate sample report.

- A multipoint calibration to confirm instrument linearity is performed for the same ten monitored hydrocarbons as in the control standard (four in the light-end HC analysis, and six in the mid-range HC analysis). It is done for new instruments, after making instrument modifications which can affect linearity, and at least once per year. The multipoint consists of at least five concentrations, each above the maximum allowable LOD (20 ppbC), about evenly distributed over the range of expected sample concentration (from about 100 ppbC to about 10000 ppbC). Each concentration is measured at least twice. A linear regression analysis is performed using concentration and average area counts to determine the regression correlation coefficient (r). The r must be greater than 0.995 to be considered sufficiently linear to ensure the accuracy of the daily one-point calibration. Figure 7 illustrates a typical multipoint calibration.
- 8.7 The LOD for the same ten monitored hydrocarbons as in the control standard must be determined for new instruments, after making modifications which can affect linearity and/or sensitivity and at least once per year. To make the calculations, it is necessary to perform a multipoint calibration consisting of at least four "low" concentration levels, each above the expected LOD. The LOD is calculated using the following equation:

$$LOD (ppbC) = (|b| + (t \times s))$$
m

where

the absolute value of the y-intercept, area counts
 the slope of the linear regression, area counts/ppbC

s = the standard deviation of at least five measurements of the lowest

concentration standard, area counts

the t-factor for 99 percent confidence for a one-sided normal (Gaussian)

distribution, dimentionless

The number of degrees of freedom is equal to the number of runs (n) of the lowest concentration standard, minus one. An abbreviated t-table is (Ref.9.1):

| Degrees of Freedom (n-1) | t-values |
|--------------------------|----------|
| 4 | 3.7 |
| 5 | 3.4 |
| 6 | 3.1 |
| 7 | 3.0 |

- 8.7.1 The concentration of the lowest standard must be greater than the calculated laboratory LOD, and not more than five times the estimated LOD. The maximum allowable LOD for each compound is 20 ppbC. The calculated laboratory LOD must be equal to or lower than the maximum allowable LOD.
- 8.7.2 All peaks identified as target compounds that are equal to or exceed the maximum allowable LOD must be reported.
- 8.7.3 If the calculated laboratory LOD is less than the maximum allowable LOD, SLB may set its reporting limit at either the maximum allowable LOD or the calculated laboratory LOD. The current reporting limit is 15 ppbC. Figure 8 shows a typical LOD determination.

9.0 **References**

- 9.1 California Environmental Protection Agency, Air Resources Board, "California Nonmethane Organic Gas Test Procedures", Part D (Method No. 1002 "Determination of C2 to C5 hydrocarbons in Automotive Source Samples by Gas Chromatography") and Part E (Method No. 1003 "Determination of C6 to C12 hydrocarbons in Automotive Source Samples by Gas Chromatography"), May 31, 1996.
- 9.2 Code of Federal Regulations, Title 40, Part 86.
- 9.3 Randy Bramston-Cook in Lotus Consulting monographs on "Speciation of Hydrocarbons in Vehicle Exhaust and Ambient Air" and "Quantitation of Hydrocarbons in Vehicle Exhaust and Ambient Air", 1996.
- 9.4 "Standard Test Method for C₁ Through C₅ Hydrocarbons in the Atmosphere by Gas Chromatography", ASTM Standards on Chromatography (1981).
- 9.5 U.S. Environmental Protection Agency, <u>Compendium of Methods for the Determination of Toxic Organic Compounds in Ambient Air (Method T93-15)</u>, EPA-600/4-89-017, Research Triangle Park, North Carolina, June 1989.

Table 2. Typical GC Operating Conditions

| GC Parameters | Light-end GC | Mid-range GC | | | | | |
|------------------------------|--------------------|-------------------|---------------------|--|--|--|--|
| | | "Fast" Analysis | "Regular" Analysis | | | | |
| Column type | Alumina/KCl PLOT | DB-1 WCOT | DB-1 WCOT | | | | |
| Column size | 0.32 mm ID, 50 m | 0.18 mm ID, 40 m | 0.32 mm ID, 60m | | | | |
| Helium carrier velocity | 30 cm/sec at 200°C | 2 mL/min at 200°C | 30 cm/sec at 200°C | | | | |
| Hydrogen gas flow(for FID) | 25 mL/min | 30 mL/min | 30 mL/min | | | | |
| "Zero" air gas flow(for FID) | 300 mL/min | 300 mL/min | 300 mL/min | | | | |
| Nitrogen make-up flow | 30 mL/min | 30 mL/min | 30 mL/min | | | | |
| Auxiliary sample valve | 150°C | 150°C | 150°C | | | | |
| Precolumn temperature | 40°C | N/A | N/A | | | | |
| Column temperature | 0°C(hold 7 min) | -50°C, hold | -50°C, (hold 2 min) | | | | |
| | 5°C/min to 50°C | 8°C/min to 170°C | 3°C/min to 170°C | | | | |
| | 10°C/min to 200°C | 50°C/min to 250°C | 47°C/min to 200°C | | | | |
| | hold 15 min | End | 30°C/min to 250°C | | | | |
| Cryotrap/Injector Program | -180°C | -180°C | -180°C | | | | |
| | 250°/min to 200°C | 250°/min to 200°C | 250°C/min to 200°C | | | | |
| | Hold for 28 min | hold for 28 min | hold for 70 min | | | | |
| Detector Temperature | 250°C | 300°C | 300°C | | | | |
| Detector Attenuation | 8 | 8 | 8 | | | | |
| Detector range | 12 | 12 | 12 | | | | |
| Total analysis time | 30 min | 35 min | 90 min | | | | |

Title : Light-End GC#4L-7992-PC307
Run File : C:\STAR\MODULE16\49C31157.RUN
Method File : C:\STAR\GC4MR.MTH Sample ID : QA C3 AAL8623 Injection Date: 10-FEB-99 8:25 AM Calculation Date: 10-FEB-99 8:42 AM Operator : LOS Detector Type: ADCB (10 Volts) Operator . LCC Workstation: AUTOGC4 Bus Address : 16 Sample Rate : 10.00 Hz Run Time : 17.112 min Instrument : MR-GC4 : A = FID LE Channel ****** Star Chromatography Workstation ****** Version 4.51 ********** Chart Speed = 1.15 cm/min Start Time = 0.000 min Attenuation = 478 Zero Offset = 3% End Time = 17.112 min Min / Tick = 1.00 1.15 cm/min Attenuation = 478 0.75 1.00 Volts 0.50 0.00 0.25 4 -+GR 5 -GR 6 +GR 9 -GR 10 11 12 +GR 300 C3 13 -GR 15 -16 <WI=4.0*

Figure 1

Title : Mid-Range GC#4M-8020-PC307 Run File : C:\STAR\MODULE16\49C31157.RUN Method File : C:\STAR\GC4MR.MTH Sample ID : QA C3 AAL8623 Injection Date: 10-FEB-99 8:25 AM Calculation Date: 10-FEB-99 8:42 AM Operator : LOS Detector Type: ADCB (10 Volts) Workstation: AUTOGC4 Bus Address : 16 Sample Rate : 10.00 Hz Run Time : 17.112 min Instrument : MR-GC4 Channel : B = FID MR******* Star Chromatography Workstation ****** Version 4.51 ********** Attenuation = 276 Zero Offset = 2% End Time = 17.112 min Min / Tick = 1.00 Chart Speed = 1.15 cm/min Attenuation = 276 0.000 min Start Time = 0.5 0.6 +GR 0.0 0.2 0.3 0.4 2.078 100 C1 2 -3 -<WI=4.0 +GR 5 -6 --GR* 300 C3 9. -GR 10 11 13 -14 15 -+GR 16 -GR 17 -

Figure 2

Title : Light-End GC#4L-7992-PC307 Run File : C:\STAR\MODULE16\49GH1106.RUN
Method File : C:\STAR\GC4MR.MTH Sample ID : ssgh Injection Date: 14-JAN-99 12:42 PM Calculation Date: 14-JAN-99 1:58 PM Operator : LOS Detector Type: ADCB (10 Volts) Bus Address : 16 Sample Rate : 10.00 Hz Workstation: AUTOGC4 Instrument : MR-GC4 Sample Rate : 76.003 min : A = FID LEChannel Run Time ******* Star Chromatography Workstation ****** Version 4.51 *********** Chart Speed = 0.73 cm/minAttenuation = 442Zero Offset = 2% End Time = 27.000 min Min / Tick = 1.00Start Time = 0.000 min 0.00 0.50 0.75 1.00 Volts 2 +GR 200 C2 4 849 5 -GR 6 7 +GR 8 250 C2= 8.596 9 -GR 10 11 12 +GR 300 C3 12.540 13 -GR 14 15 350 C3= -15.805 16 <WI=4.0* 395 C2yne 400 C4 18 19 455 1C4= -19.576 19.864 20 21 500 C5 -21.408 22 23 24 <WI=2.0 600 C6 -24.883 +11 26

Figure 3

27

Figure 4

ETHANE




Figure 5

REPLICATES ANALYSIS

| _ | | STATUS | ¥ | A | Ν | ¥ | ¥ | ¥ | FAIL | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ |
|----------|-------------|-----------|----------|-----------------------|------------|------------|-------------|------------|------------|----------------|------------|------------|------------|------------|------------|------------|-----------------|------------|------------|------------|------------|
| 111 | | MAX. | 2 | 15 | 5 | 20 | 5 | 8 | 8 | 용 | 6 | 용 | 용 | 6 | 용 | R | 5 | 8 | 6 | 6 | 5 |
| PROPANE | 200 | %RPD | Б | 2 | 4 | = | വ | - | R | Ξ | | Ξ | 2 | 2 | 7 | 2 | 14 | 2 | D. | Ξ | 6 |
| <u>.</u> | <i>=001</i> | RUN #2 | £ 52 | 134531 | 279 | 125 | 3 91 | 9 | 8 | ය | 24 | 8 | ස | 43 | 2 | 8 | 538 | 222 | 44 | æ | 583 |
| | | BUN # | 38 | 137927 | 230 | 140 | 248 | ē | 113 | න | 92 | 72 | 25 | 42 | 8 | 8 | 23 3 | 526 | 42 | 34 | 315 |
| _ | | STATUS | ¥ | ŏ | ¥ | ă | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ | ¥ |
| | | MAX. | 2 | ##### | 15 | 15 | 5 | 15 | 15 | 5 | 6 | 8 | 8 | 15 | 5 | 5 | 5 | 5 | 5 | 5 | 15 |
| ETHENE | 200 | %RPD | 52 | #VALUE! | 0 | 0 | 0 | - | 2 | 2 | 0 | 00 | - | - | - | - | 0 | - | - | 0 | വ |
| | <i>=001</i> | RUN #2 | ‡ ‡ | Ą | 648 | 10849 | 22650 | 809 | 1380 | 573 | 13 | ន | 127 | 1613 | 283 | 222 | 829 | 1830 | 743 | 623 | 14703 |
| | | BUN #1 | 9 9 | Ą | 650 | 10842 | 22611 | 604 | 1410 | 261 | 13 | ਹ | 128 | 1623 | 227 | 559 | 929 | 1845 | 738 | 622 | 15390 |
| - | | TEST # | 2C24 BG2 | ROPANE CYL | 10C24 BG1 | 235C1 bag2 | 235U1 bag2 | 10C26 bag1 | 13C25 bag1 | 10C27 bag1 | 10C28 bag3 | 13C28 bag3 | 13C29 bag3 | 13C30 bag1 | 10C32 bag1 | 10C33 bag1 | 6C1 bag1 | 5C1 bag1 | 6C3 bag1 | 6C4 bag1 | 262C1 BG2 |
| | | DATE | 8 | 10/05/1998 PROPANE CY | 10/15/1998 | 10/16/1998 | 10/19/1998 | 10/20/1998 | 10/22/1998 | 10/23/1998 | 10/27/1998 | 10/28/1998 | 10/29/1998 | 10/30/1998 | 11/03/1998 | 11/04/1998 | 11/05/1998 | 11/06/1998 | 11/10/1998 | 11/20/1998 | 12/16/1998 |

Figure 6

| | SUMARY | RESULTS OF LINE | ARITY CHECK | | |
|---|-------------------------|-----------------|------------------------------|----------------------------|------------------------|
| Period: Gas Chromatograph: | 03/25/1996 No. 6 | (MR=8018) | | | |
| HYDROCARBON: | HEXANE | | | | R= 0.9978 M= 138.50 |
| Standard ID | Dilution Factor | Number of Runs | Average Area Count | ppmC W/R to List Values | 141- 150.50 |
| NIST STDAAL10318(v v) | Undiluted | 5 | 42168 | 274 | |
| NIST STDAAL111183(w w) NIST STDAAL111183(w w) | 0.1 0.2 | 3 3 | 140478 305748 | 888 1776 | |
| NIST STDAAL111183(w w) | 0.4 | 3 | 544257 | 3552 | |
| NIST STDAAL111183(w w) NIST STDAAL111183(w w) | 0.7 0.9 | 3 3 | 914379 1170017 | 6216 7992 | |
| NIST STDAAL111183(w w) | Undiluted | 6 | 1211005 | 8880 | |
| HYDROCARBON: | BENEZENE | | | | R= 0.9987 |
| Standard ID | Dilution Factor | Number of Runs | Average Area | ppmC W/R to | M= 140.74 |
| NIST STDAAL10318(v v) | Undiluted | 5 | Count 41058 | List Values 258 | |
| NIST STDAAL111183(w´w) | 0.1 | 3 | 146957 | 1002 | |
| NIST STDAAL111183(w w) | 0.2 | 3 | 342511 | 2004 | |
| NIST STDAAL111183(w w) NIST STDAAL111183(w w) | 0.4 0.7 | 3 3 | 594379 1025474 | 4008 7014 | |
| NIST STD.:AAL111183(w w) | 0.9 | 3 | 1322584 | 9018 | |
| NIST STDAAL111183(w w) | Undiluted | 6 | 1394944 | 10020 | |
| HYDROCARBON: | TOLUENE | | | | R= 0.9563 M= 157.96 |
| Standard ID | Dilution Factor | Number of Runs | Average Area Count | ppmC W/R to List Values | 101.50 |
| NIST STDAAL10318(v v) | Undiluted | 5 | 39888 | 218 | |
| NIST STDAAL111183(w w) | 0.1 0.2 | 3 3 | 153211 376368 | 1008 2016 | |
| NIST STDAAL111183(w w) NIST STDAAL111183(w w) | 0.4 | 3 | 627248 | 4032 | |
| NIST STDAAL111183(w w) | 0.7 | 3 | 1085002 | 7058 | |
| NIST STDAAL111183(w w) NIST STDAAL111183(w w) | 0.9 Undiluted | 3 6 | 1389105 1500967 | 6072 10080 | |
| HYDROCARBON: | OCTANE | | | | R= 0.9990 |
| Standard ID | Dilution Factor | Number of Runs | Average Area | ppmC W/R to | M= 149.20 |
| | | | Count | List Values | |
| NIST STDAAL10318(v v) | Undiluted | 5 | 40931 | 278 | |
| NIST STDAAL111183(w w) NIST STDAAL111183(w w) | 0.1 0.2 | 3 3 | 153096 355672 | 1000 2000 | |
| NIST STDAAL111183(w w) | 0.4 | 3 | 621510 | 4000 | |
| NIST STDAAL111183(w w) | 0.7 | 3 | 1087118 | 7000 | |
| NIST STDAAL111183(w w) NIST STDAAL111183(w w) | 0.9 Undiluted | 3 6 | 1381968 1477929 | 9000 10000 | |
| | | | | .5566 | D 6 666 |
| HYDROCARBON: | XYLENE | | | C.W.D. | R= 0.9986 M= 151.05 |
| Standard ID | | Number of Runs | Average Area Count | ppmC W/R to List Values | |
| NIST STDAAL10318(v v) NIST STDAAL111183(w w) | Undiluted 0.1 | 5 3 | 39835 149961 | 253 992 | |
| NIST STD.:AALTTT163(w w) | 0.1 | 3 | 384404 | 1984 | |
| NIST STDAAL111183(w w) | 0.4 | 3 | 596373 | 3968 | |
| NIST STDAAL111183(w w) NIST STDAAL111183(w w) | 0.7 0.9 | 3 3 | 1058690 1381909 | 6944 8928 | |
| NIST STD.: AAL111183(w w) | Undiluted | 6 | 1512694 | 9920 | |
| HYDROCARBON: | DECANE | | | | R= 0.9989 M= 146.15 |
| Standard ID | Dilution Factor | Number of Runs | Average Area | ppmC W/R to | W- 140.13 |
| NIST STDAAL10318(v v) | Undiluted | 5 | Count 39689 | List Values 242 | |
| NIST STDAAL111183(w´w) | 0.1 | 3 | 138756 | 936 | |
| NIST STDAAL111183(w w) | 0.2 | 3 | 338353 | 1872 | |
| NIST STDAAL111183(w w) = | 0.4 | 3 | 540051 | 3744 | |
| | | | 990061 | | |
| NIST STDAAL111183(w w) NIST STDAAL111183(w w) NIST STDAAL111183(w w) NIST STDAAL111183(w w) | 0.7 0.9 Undiluted | 3 3 6 | 990061 1266139 1370613 | 6662 8424 9360 | |

Figure 7

| | SUMA | RY RESULTS OF L | OD CHECK | | |
|--|---------------------|-----------------|-----------------------|-------------------|------------------------|
| Period: Gas Chromatograph: | 03/25/1996 No. 6 | (MR=8018) | | | |
| HYDROCARBON: | HEXANE | | LOD= | 5.77 ppbc | R= 0.9884 |
| Standard ID | Dilution Factor | Number of Runs | Average Area | ppmC W/R to | M= 145.51 |
| NUCT CTD 0 01 102104 | 0.4 | F | Count | List Values | |
| NIST STDAAL10318(v v) NIST STDAAL10318(v v) | 0.1 0.2 | 5 3 | 4396 8574 | 27 55 | |
| VIST STD.: AAL10316(V V) | 0.5 | 3 | 21375 | 137 | |
| VIST STDAAL10318(v v) | 0.8 | 3 | 28397 | 219 | |
| NIST STDAAL10318(v v) | Undiluted | 3 | 42736 | 274 | |
| HYDROCARBON: | BENZENE | | LOD= | 2.7 ppbc | R= 0.9847 |
| Standard ID | Dilution Factor | Number of Runs | Average Area | ppmC W/R to | M= 151.32 |
| | | | Count | List Values | |
| NIST STDAAL10318(v v) | 0.1 | 5 | 4247 | 26 | |
| NIST STDAAL10318(v v) | 0.2 | 3 | 8073 | 52 | |
| NIST STDAAL10318(v v) | 0.5 | 3 | 20407 | 129 | |
| NIST STDAAL10318(v v) | 0.8 | 3 | 26716 | 205 | |
| NIST STDAAL10318(v v) | Undiluted | 3 | 41950 | 258 | |
| HYDROCARBON: | TOLUENE | | LOD= | 1.9 ppbc | R= 0.9876 |
| Standard ID | Dilution Factor | Number of Runs | Average Area | ppmC W/R to | M= 188.33 |
| | | | Count | List Values | |
| NIST STDAAL10318(v v) | 0.1 | 5 | 41.62 | 22 | |
| NIST STDAAL10318(v v) | 0.2 | 3 | 7646 | 44 | |
| NIST STDAAL10318(v v) | 0.5 | 3 | 19789 | 109 | |
| NIST STDAAL10318(v v) | 0.8 | 3 | 26792 | 174 | |
| NIST STDAAL10318(v v) | Undiluted | 3 | 40478 | 218 | |
| HYDROCARBON: | OCTANE | | LOD= | 2.21 ppbc | R= 0.9895 M= 139.87 |
| Standard ID | Dilution Factor | Number of Runs | Average Area | ppmC W/R to | 155.61 |
| | | | Count | List Values | |
| NIST STDAAL10318(v v) | 0.1 | 5 | 4254 | 28 | |
| NIST STDAAL10318(v v) | 0.2 | 3 | 8144 | 56 | |
| NIST STDAAL10318(v v) | 0.5 | 3 | 20143 | 139 | |
| NIST STDAAL10318(v v) | 0.8 | 3 | 27617 | 222 | |
| NIST STDAAL10318(v v) | Undiluted | 3 | 41409 | 278 | |
| HYDROCARBON: | XYLENE | | LOD= | 2.73 ppbc | R= 0.9829 M= 148.02 |
| Standard ID | Dilution Factor | Number of Runs | Average Area Count | ppmC W/R to | 170.02 |
| NIST STDAAL10318(v v) | 0.1 | 5 | 4295 | List Values 25 | |
| NIST STDAAL10318(v v) | 0.1 | 3 | 7579 | 25 51 | |
| NIST STD.: AAL10316(v v) | 0.5 | 3 | 19235 | 126 | |
| NIST STDAAL10318(v v) | 0.8 | 3 | 25572 | 202 | |
| NIST STDAAL10318(v v) | Undiluted | 3 | 40529 | 253 | |
| HYDROCARBON: | DECANE | | LOD= | 4.43 ppbc | R= 0.9846 |
| Standard ID | Dilution Factor | Number of Runs | Average Area | ppmC W/R to | M= 156.06 |
| | | | Count | List Values | |
| NIST STDAAL10318(v v) | 0.1 | 5 | 4620 | 24 | |
| NIST STDAAL10318(v v) | 0.2 | 3 | 8282 | 48 | |
| NIST STDAAL10318(v v) | 0.5 | 3 | 18427 | 121 | |
| NIST STDAAL10318(v v) | 0.8 | 3 | 26620 | 194 | |
| NIST STDAAL10318(v v) | Undiluted | 3 | 41190 | 242 | |

Figure 8