

Novel Approaches for Speciation of Platinum and Vanadium in Mobile Source Emissions


**Martin Shafer¹, James Schauer¹, Water Copan²,
Alberto Ayala³, Shaohua Hu³, Jorn Herner³**

**¹University of Wisconsin-Madison
Environmental Chemistry & Technology**

²Clean Diesel Technologies Inc.

³California Air Resources Board

AAAR 27TH Annual Conference



Motivation

- Controlling emissions from mobile sources are critical for continued reduction in health impacts of air pollution, and for addressing regional and global climate impacts.
- Most current and proposed emission control strategies for diesel and gasoline engines employ metal catalysts to reduce tailpipe emissions of regulated species.
 - ❖ *Gasoline Three-Way-Catalysts (Pt, Pd, Rh)*
 - ❖ *Diesel Fuel Based Catalysts (Pt-FBC)*
 - ❖ *Diesel Particulate Filters (Pt-Catalyzed)*
 - ❖ *Diesel Selective Catalytic Reactors (V-SCR)*
- The use of these metals raises concerns about potential environmental contamination and the health implications of widespread trace metal dissemination.

MOTIVATION

- The toxicological responses of many metals (e.g. Cr, Cu, Mn, Pt, V) are determined by the specific chemical & physical speciation in the emissions.
 - ❖ **Platinum**: soluble, oxidized, halogenated (e.g. chloroplatinic acids) species are 500 fold more toxic than metallic species.
 - ❖ **Vanadium**: pentoxide V(V) species exhibits much greater toxicity than the lower oxidation state species.
- Problem and Challenge:
 - ☹ **Extant modern methodologies provide little relevant speciation information.**
 - ☹ **Traditional techniques that are speciation capable lack the required sensitivity, particularly in the context of lower emissions from vehicles equipped with modern control devices.**

Engine Exhaust PM Characterization

□ Elemental & Isotopic Characterization

- *Magnetic sector (high resolution) ICP-MS*
- *Sensitivity and Interference Isolation*

□ Chemical & Physical Speciation of PM

- Solubility
- Oxidation State

▪ Soluble Species:

- *Long-path (100 cm) Spectrophotometry with Characteristic Ligand e.g. Fe(II)/Fe(III), Mn(II)/Mn(>III), Cr(III)/Cr(VI)*
- *Immobilized Ligand – Selective Extraction and Elution*

▪ Direct Solids:

- *Synchrotron X-ray Absorption Spectroscopy*
- Colloid Charge: *Ion Chromatography (DEAE, SAX micro-columns) → ICP-MS*
- Colloid Size: *Ultrafiltration (10, 100 kDa) → ICP-MS*

Complementary Total and Extractable Methods

Speciation Background

- **Platinum**

- Oxidation States (0, II, IV). Group 8 transition metal.
- Higher oxidation states may be more soluble.
- Chloroplatinic/um salts (H, NH₄, K, Na) are very soluble.

- **Vanadium**


- [Ar] 3d³ 4s². Forms oxyanions. Amphoteric. Redox Active.
- Oxidation States: (0,II,III,**IV,V**).
- Higher oxidation states more soluble in water due to hydrolysis.
- **V(V)**: (high pH) VO₄³⁻, HVO₄²⁻, H₂VO₄⁻, H₃VO₄, VO₂⁺ (low pH).
- **V(IV)**: cationic (VO²⁺)

Platinum: Sources and Receptors Under Study

- ❑ Three-Way-Catalysts
- ❑ Size-Resolved PM from Engines burning Platinum-Amended Diesel Fuel
- ❑ PM from Platinum-Catalyzed DPF
- ❑ Roadside Dust / Soils




25 mm PCIS Substrate



**Gasoline
Engine
Catalytic
Converter**

What is a Fuel-Borne Catalyst?

- **Catalyst dosed directly into diesel fuel**
 - Pt / Ce fuel-soluble bimetallic catalyst
 - delivered *in situ*
- **Active in high temperature combustion zone**
 - higher efficiency of fuel HC combustion
- **FBC intimate contact with PM**
 - more complete combustion of solid C, HC
 - uniformly distributed across PM size range
 - no increase in ultra-fines
- **Delivers Catalyst to DOC / DPF**
 - fresh catalyst surface replenishment
 - same active forms
 - permits lower lifetime use of Pt



Platinum: Analytical Speciation of Engine PM

■ Solid Phase

- Total: microwave/acid digestion – HR-ICP-MS
- **Oxidation State: Synchrotron XAS**

■ Extractable Species

- Primarily oxidized and halogenated species
- Total “soluble” - HR-ICP-MS
 - ❖ *Water*
 - ❖ *Methanol*
 - ❖ *DCM*
- Ultra-filtration: colloidal versus “dissolved”
- Ion Chromatography: anionic versus cationic
- HPLC-IC-HR-ICP-MS: specific chemical species

■ Particle Size Distribution

- Sioutas PCIS (5 size-cuts)

Analytical Challenge: <1 ng extractable platinum (10-30 pg in specific fractions).

Speciated Water Soluble Platinum in Diesel PM


FTP Cycle Means

Extractable fraction = < 3%.

Large colloidal fraction (44% of extractable species).

Dissolved (<10 kDa) species exhibit significant anionic character on DEAE (42%).

Shafer, M.M., J.J. Schauer, W. Copan, J. Peter-Hoblyn, B. Sprague, and J. Valentine. 2007. Investigation of platinum and cerium from use of a fuel-based catalyst. *SAE 2006 Transactions Journal of Fuels and Lubricants* – 2006-1-1517:491-503.




Vanadium: Analytical Speciation of Engine PM and Urban Aerosol


- Solid Phase
 - Total: microwave/acid digestion – HR-ICP-MS
 - **Oxidation State: Synchrotron XAS**
- Extractable Species
 - Primarily V(V) and V(IV)
 - Total water and acetate soluble: HR-ICP-MS
 - **[V(II), V(III)] V(IV), V(V) : Immobilized ligand speciation**
- Particle Size Distribution
 - Sioutas PCIS (5 size-cuts)

Significant Analytical Challenge: 0.2-2 ng total vanadium from Dyno Trials for speciation studies.


Vanadium Catalyst: Total Vanadium




Vanadium Catalyst: Soluble Vanadium



SCRT® System Schematic



Vanadium Catalyst: Soluble V (% of Total)



Vanadium Oxidation State Speciation with Chelex

Chelex: cation exchanger at higher pH, anion exchanger at lower pH.

In pH range of 4 to 7.4 both cations and anions are adsorbed.

■ Preparation & Extraction

- ❖ Micro-columns of Chelex (immobilized iminodiacetate)
 - 0.2 g of perchloric acid cleaned and acetate buffered (pH=4.5) resin
- ❖ Samples extracted in 2 mM sodium acetate buffer
 - 1.5 mL nitrogen purged buffer in purged cryo-vial
 - 60 minutes with defined agitation (under nitrogen canopy)


■ Separation – Speciation

- ❖ 1.0 mL of sample loaded onto column
- ❖ Process through column at 1.0 mL per minute – collect fraction
- ❖ Elute column with 4 x 1 mL of 0.1 M ammonium hydroxide – collect fraction = **V(V)**.
- ❖ Elute column with 2 x 1 mL of 0.2 M perchloric acid – collect fraction = **V(IV)**

■ Vanadium Quantification

- ❖ Magnetic Sector ICP-MS in medium resolution with on-line standard addition
- ❖ 10,000 cps/ppb V. Background = ~2 cps (0.2 ng L⁻¹). 1-5% RSD.

Vanadium Oxidation State Speciation with Chelex




0.044 ± 0.037 ng

0.016 ± 0.011 ng


(IV) DL=0.11 ng
(0.44 ppm)

(V) DL=0.03 ng
(0.14 ppm)


Vanadium Speciation – Defined Standards



Vanadium Speciation – Environmental Matrices




Vanadium Speciation – Engine Exhaust PM



ADVANCES



1. Micro-scale
2. DLs improved by 10-50x
3. Coupling to HR-ICP-MS



Wang D. and S. Sanudo-Wilhelmy 2008 Marine Chemistry.


Synchrotron X-Ray Absorption Spectroscopy

- Direct Solids Analysis – complementary to solution phase tools.
- XANES – (oxidation state) and/or EXAFS – (nearest neighbor chemical bonding environment). XRD – (characteristic diffraction)
- Range of spatial scales (with/without micro-focused beamline).



Majestic, B.J., J.J. Schauer, M.M. Shafer. 2007. Application of Synchrotron Radiation for Measurement of Iron Red-Ox Speciation in Atmospherically Processed Aerosols. *Atmospheric Chemistry and Physics* 7:2475-2487.


EXAFS Spectra of Platinum Reference Materials




LBL-Advanced Light Source




Pt-foil = platinum foil; Pt-alumina = 5% platinum on alumina; Pt-C = 5% platinum on graphite; Pt(II)Cl₂ = platinum(II) chloride; TCP = potassium platinum(II) tetrachloroplatinate; PtO₂ = platinum(IV) oxide; Pt(IV)Cl₄ = platinum(IV) chloride; HCP = potassium platinum(IV) hexachloroplatinate




A



C



D



B

Pt speciation was studied in a 4 year old 3-way automobile catalyst.

A 30 μm thick, quartz slide mounted, longitudinal section of the center of the catalyst was prepared.

(A) Light microscope image; outlined area was examined with XRF mapping.

(B) Red (Pt) – green (Cu) – blue (Ce) XRF-derived tricolor map. Pt L3-edge extended-XANES spectra were collected at spots 0-2.


(C) The e-XANES spectra (11,466-12,077 eV) were fit with reference spectra - Pt foil, 5% Pt in alumina matrix, 5 % Pt in carbon matrix, Pt(II)Cl_2 , Pt(IV)Cl_4 , PtO_2 , $\text{K}_2\text{Pt(IV)Cl}_6 \cdot \text{H}_2\text{O}$, and $\text{K}_2\text{Pt(II)Cl}_4$ – by linear least squares method

(D) Select reference spectra and an example fit shown in C and D).

Modeled fraction of oxidized Pt is significant.

XRF Map, and Extended-XANES Fitting Results, of Diesel Exhaust Particulates Trapped on a Diesel Particulate Filter (engine running a Pt-FBC)

X-ray Fluorescence Map




Spot #	Elemental Composition	Linear Least-Squares Fit Results
0	Ce, Zn, S	88% Pt-alumina, 7% Pt-foil, 5% TCP
1	Ce	56% Pt-alumina, 21% Pt-C, 24% PtO ₂
2	Fe-rich	87% Pt-alumina, 10% PtO ₂ , 3% HCP
3	Pt alone	81% Pt-alumina, 18% PtO ₂ , 2% TCP
4	Ca, S	59% Pt-foil, 41% PtO ₂
5	S-rich	86% Pt-C, 11% TCP, 7% Pt-alumina
6	Ce	86% Pt-alumina, 14% PtO ₂ , 3% HCP
7	Pt-rich	100% Pt-foil
8	Ce, Pt-rich	64% Pt-alumina, 11% PtO ₂ , 26% Pt-foil
9	Ce-rich	36% Pt-foil, 7% HCP, 59% Pt-C

Large heterogeneities in particle composition are observed with many particles exhibiting a significant oxidized platinum component.

Strong evidence for PtO₂ (14-25% in many spots, up to 40% when associated with Ca and S), and less compelling evidence for presence of chloroplatinates.

Platinum bulk-XANES

XAS Spectra of Platinum Standards and Diesel Engine Exhaust PM (Running Pt-FBC Fuel)



ARGONNE
NATIONAL LABORATORY




Shafer, M.M., J.J. Schauer, W. Copan, J. Peter-Hoblyn, B. Sprague, and J. Valentine. 2007. Investigation of platinum and cerium from use of a fuel-based catalyst. *SAE 2006 Transactions Journal of Fuels and Lubricants* – 2006-1-1517:491-503.



University of Wisconsin - Madison
Wisconsin State Laboratory of Hygiene

XRF Maps of Diesel Exhaust Particulates Impacted on Teflon PCIS Substrates (engine running with a **vanadium SCR**)



Strong evidence of V, Ti –rich particles in engine PM

Preliminary e-XANES suggests V(V)

Acknowledgments



HEI



Comparison of Water & Solvent – Speciated Platinum Emissions

- ❑ Significantly more Pt extracted with MeOH ($28 \pm 1.4\%$) than with DCM ($0.79 \pm 0.15\%$) or water.
- ❑ DCM should not extract ionic Pt species (except via physical process).
- ❑ DCM more selective in isolating any Pt associated with organic matter.
- ❑ Methanol (MeOH) extracts and disperses both polar and non-polar species. Breaks up diesel PM soot matrix. High MeOH extractables fraction does not indicate the presence of a large pool of organo-Pt species.

